

Oxygen Concentrators

Perfecto2 (IRC5P, IRC5PO2)

Perfecto2V (IRC5PO2V)

PerfectoW (IRC5PW, IRC5PO2W)

Platinum 5 (IRC5LX)

Non Platinum 5LX

Platinum 10 (IRC10LX, IRC10LXO2)

HomeFill™ Oxygen System

DODD

Check 02 Plus Conc. Analyzer



Yes, you can.

Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 93-38. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-333-6900

For Faxed In Orders Call:

1-800-678-4682 or

1-440-366-0709

For Technical Service Call:

1-800-832-4707

Literature Can Be Found at:

www.invacare.com or

1-800-828-6272

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers

© 2020 Invacare Corporation, Elyria, Ohio. Invacare, the Medallion Design, Venture and Yes, you can. are registered trademarks of the Invacare Corporation. Platinum, Check O2 Plus and HomeFill are trademarks of Invacare Corporation. Palnut is a registered trademark of Transtechnology Corporation. Tygon is a registered trademark of Norton Company. Information is subject to change without notification.

TABLE of CONTENTS

Description	Page Number
Accessories - Concentrator Tool Kit	92
Accessories - Venture Demand Oxygen Delivery Device & Cylinder	90
Cabinet and Base (Perfecto2)	16
Cabinet and Base (Perfecto2V)	18
Cabinet and Base (PerfectoW)	20
Cabinet and Base (Platinum 10)	14
Cabinet and Base (Platinum 5 Models)	12
Cabinet, Base, Cylinder and Tabletop (Non Platinum 5 Models)	9
Compressor and Mounting Hardware (Non Platinum 5 Models)	22
Compressor and Mounting Hardware (Perfecto2/PerfectoW)	31
Compressor and Mounting Hardware (Perfecto2V)	33
Compressor and Mounting Hardware (Platinum 10 Models)	28
Compressor and Mounting Hardware (Platinum 5 Models)	25
Control Panel and P.C. Board (Platinum 5 Models - After 12/5/02)	44
Control Panel and P.C. Board (Platinum 5 Models - Before 12/5/02)	41
Control Panel and P.C. Board (Non Platinum 5 Models - After 12/14/00)	38
Control Panel and P.C. Board (Non Platinum 5 Models - Before 12/15/00)	35
Control Panel and P.C. Board (Perfecto2, Perfecto2V, PerfectoW)	50
Control Panel and P.C. Board (Platinum 10)	47
HomeFill Base and HomeFill Stand Components	86
HomeFill Cylinder Assembly	88

TABLE of CONTENTS

Description	Page Number
Invacare Check 02 Plus Concentrator Analyzer	85
Miscellaneous Components	94
Packaging, Literature, and Preventative Maintenance Kit	4
Quick Reference Chart	6
Remote Flowmeters	82
Sieve Bed Assembly and Product Tank (Non Platinum 5 Models)	66
Sieve Bed Assembly and Product Tank (Perfecto2)	74
Sieve Bed Assembly and Product Tank (Perfecto2V)	77
Sieve Bed Assembly and Product Tank (PerfectoW)	79
Sieve Bed Assembly and Product Tank (Platinum 10)	72
Sieve Bed Assembly and Product Tank (Platinum 5 Models)	69
Sound Box Assembly (Non Platinum 5 Models)	52
Sound Box Assembly (Perfecto2, Perfecto2V)	62
Sound Box Assembly (PerfectoW)	64
Sound Box Assembly (Platinum 10)	59
Sound Box Assembly (Platinum 5 Models)	56
Warranty Replacement Parts	95

Packaging, Literature, and Maintenance Kits

Literature and Preventative Maintenance Kit

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A,B	1106666-NLA	NLA - Owner's Manual, Domestic U.S. (IRC3LX,5LX,6LX)	1	1
	C,A	1110514-NLA	NLA - Owner's Manual, Domestic U.S. (IRC3LX,5LX)	1	1
	A,B	1059390-NLA	NLA - Owner's Manual, Canadian (IRC3LX,5LX,6LX)	1	1
	F,A	1193323	Owner's Manual, 10LXO2	1	1
	G,A	1118353-NLA	NLA - Owner's Manual, 10L	1	1
	A	1100873	Owner's Manual, Domestic U.S., HomeFill Oxygen System	1	1
	A	1114802-NLA	NLA - Owner's Manual, Canadian, HomeFill Oxygen System	1	1
	A	1075286-NLA	NLA - Operator's Manual, IOH100	1	1
	A,B	2000895	Kit, Preventative Maintenance (IRC3LX,5LX)	1	1
	D,A	1117149	Kit, Receptacle Fitting Brass	1	1
			NLA - Fitting, Brass (1/4")	1	1
			Instruction Sheet	1	1
	E,A	1123520	Kit, Preventative Maintenance, Platinum 5	1	1
		See Chart	Packaging and Literature	1	1
	F,A	1193322	Owner's Manual, Perfecto2V	1	1
	G,A	1143482-NLA	NLA - Owner's Manual, Perfecto2/Perfecto2V/PerfectoW	1	1
	A	1154302-NLA	NLA - Placard, HomeFill Compressor Guide For Facility Use	1	1
<p>NOTE: A - Not Shown. B - All Models Except Platinum 5. C - Platinum 5 Models. D - Kit includes: (1) each #1105572 Receptacle Fitting, (1) each #1117150 Instruction Sheet. E - Kit includes: (2) each #1107412 Cabinet Filter, (1) each #1131249 Inlet Filter, (1) each #P0750X030 HEPA Filter and (1) each #1134864 Preventative Maint. Record Sheet. F - Use for Concentrators built (After 9/30/16). G - Use for Concentrators built (Before 10/01/16).</p> <p>NOTE: IRCSA - Extended Warranty for IRC3LX Models. NOTE: When adjusting or replacing parts, refer to your Owner's Manual for assistance.</p>					

Packaging, Literature, and Maintenance Kits

Packaging

	Packaging Kit	Carton	Kit Top/Bottom	Top Insert	Poly Bag	Bottom Insert	Misc. Pieces Assembly
All Models Except Platinum 5	2002771 See Note(s) A	2000725-NLA	N/A	1112247-NLA	2000454	1112246-NLA	1110580-NLA See Note(s) F
Platinum 5 and 10L Models	1133694 See Note(s) H	60001254 See Note(s) J	1135509 See Note(s) D	1116624-NLA See Note(s) C	2000454	1116625-NLA See Note(s) C	N/A
10 L Model	1110846	60001254 See Note(s) J	1135509	N/A	2000454	N/A	N/A
Homefill	1102315 See Note(s) B	6000252 See Note(s) J	N/A	N/A	N/A	N/A	N/A
Perfecto2, Perfecto2V, PerfectoW	1161123 See Note(s) I	1194433 See Note(s) J	N/A	1162044	N/A	N/A	N/A
PL6 Homefill Cylinder	N/A	1151765	1151765 See Note(s) E	N/A	N/A	N/A	N/A
P9 Homefill Cylinder	N/A	1151766	1151766 See Note(s) E	N/A	N/A	N/A	N/A
PL4 Homefill Cylinder	N/A	1151767	1151767 See Note(s) E	N/A	N/A	N/A	N/A

NOTE: A - Includes Carton, Inserts, Bag, Owner's Manual, Tygon Tubing and Humidifier Adaptor.
 B - Includes Shipper Carton, Bag and Foam Inserts.
 C - Individual item no longer available, refer to Kit Part #1135509 as a replacement.
 D - Kit includes Top and Bottom Foam Inserts.
 E - Kit includes Carton and Foam Inserts.
 F - Includes Humidifier Bottle Adapter, Owner's Manual, Humidifier Bottle Tube, and Poly Bag.
 H - Kit includes Carton, Top Insert, Bottom Insert and Poly Bag
 I - Kit includes Carton, Lower Insert, and Tray.
 J - Tech Only Carton.

Quick Reference Chart

When ordering 1130666, you must order 1186230 also.

If unit was built before 01/01, then the entire Manifold/Valve Assembly Kit Part #1105552 is required for Non Platinum Concentrators. See category titled "Sound Box Assembly."

Model	Compressor Assembly	4-Way Valve	P.C. Board	Sieve Beds
3LX	1101715	1101141 See Note(s) G	1078183-NLA See Note(s) A	1101707 See Note(s) B, O
5LX Before 04/02	1101715	1101141 See Note(s) G	1078183-NLA See Note(s) A	1101707 See Note(s) B, O
5LX Platinum	1106690 See Note(s) C	1101141 See Note(s) G	1078183-NLA See Note(s) A, F	1110783-NLA See Note(s) B
5LX Platinum	1106690E-NLA See Note(s) C	1101141 See Note(s) G	1078183-NLA See Note(s) A, F	1110783-NLA See Note(s) B
5LX Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108720-NLA See Note(s) E	1110783-NLA See Note(s) B
5LX Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108720E-NLA See Note(s) E	1110783-NLA See Note(s) B
5LX02 Before 7/98	1101715	1101141 See Note(s) G	1078181-NLA See Note(s) A	1101707 See Note(s) B, O
5LX02 7/98 - 03/02	1101715	1101141 See Note(s) G	1113693 See Note(s) A	1101707 See Note(s) B, O
5LX02 Platinum	1106690 See Note(s) C	1101141 See Note(s) G	1113693 See Note(s) A, F	1161699 See Note(s) B, O
5LX02 Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108731 See Note(s) A	1161699 See Note(s) B, O
5LX02 Platinum	1157138 See Note(s) D	1101141 See Note(s) G	1108731E-NLA See Note(s) A	1161699 See Note(s) B, O
5SLX02T	1101715	1101141 See Note(s) G	1079444-NLA See Note(s) A	1101707 See Note(s) B, O
6LX02 Before 7/98	1075071-NLA	1101141 See Note(s) G	1078181-NLA See Note(s) A	1101707 See Note(s) B, O
6LX02 After 6/98	1075071-NLA	1101141 See Note(s) G	1113693 See Note(s) A	1101707 See Note(s) B, O
6LX02T	1075071-NLA	1101141 See Note(s) G	1079444-NLA See Note(s) A	1101707 See Note(s) B, O
6LX02H	1075071-NLA	1101141 See Note(s) G	1078181-NLA See Note(s) A	1101707 See Note(s) B, O
10LX	1130666	1101141 See Note(s) G	1128737	1144034 See Note(s) B
10LX	1130666E See Note(s) L	1101141 See Note(s) G	1128737	1144034 See Note(s) B
10LX02	1130666	1101141 See Note(s) G	1128742 See Note(s) M	1144034 See Note(s) B
10LX02	1130666	1101141	1194797 See Note(s) M	1144034
10LX02	1130666E See Note(s) L	1101141 See Note(s) G	1128742 See Note(s) M	1144034 See Note(s) B
Perfecto2 (IRC5P), PerfectoW (IRCPW) w/o SensO2	1162038 See Note(s) H	1101141	1143475	1161257 See Note(s) I

Quick Reference Chart

When ordering 1130666, you must order 1186230 also.

If unit was built before 01/01, then the entire Manifold/Valve Assembly Kit Part #1105552 is required for Non Platinum Concentrators. See category titled "Sound Box Assembly."

Model	Compressor Assembly	4-Way Valve	P.C. Board	Sieve Beds
Perfecto2 (IRC5PO2), PerfectoW (IRCPO2W) w/ SensO2	1162038 See Note(s) H	1101141	1143486	1161257 See Note(s) I
Perfecto2 (IRC5PO2), PerfectoW (IRCPO2W) w/ SensO2	1162038 See Note(s) H	1101141	1143486E	1161257 See Note(s) I
Perfecto2V (IRC5PO2V) w/ SensO2 See Note(s) O	1165097 See Note(s) J	1101141	1164808-NLA See Note(s) N	1165099 See Note(s) K
Perfecto2V (IRC5PO2V) w/ SensO2	1165097 See Note(s) J	1101141	1194768 See Note(s) N	1165099 See Note(s) O, K
<p>NOTE:</p> <p>A - Includes Instruction Sheet.</p> <p>B - Includes (2) Sieve Beds, (2) Clamps, (2) Hose Clamps and (4) Tie Wraps.</p> <p>C - Before 4/14/04.</p> <p>D - After 4/13/04.</p> <p>E - After 11/02. Part #1108720 is NLA. See part #1108720E as a replacement.</p> <p>F - Before 12/02.</p> <p>G - If unit was built before 01/01, then the entire Manifold/Valve Assembly Kit Part #1105552 is required for Non Platinum Concentrators. See category titled "Sound Box Assembly."</p> <p>H - When replacing the compressor assembly, part number 1131250 Capacitor w/ Tie Wraps is required also.</p> <p>I - Includes (2) Sieve Beds, (2) Tie Wraps 8" Long, (4) Tie Wraps 17-1/2" Long</p> <p>J - When replacing the compressor assembly, part number 1165098 Capacitor w/ Tie Wraps is required also.</p> <p>K - Includes (2) Sieve Beds, (4) Tie Wraps 8" Long, (4) Tie Wraps 17-1/2" Long</p> <p>L - When ordering 1130666E, you must order Kit part number 1186230.</p> <p>M - 1128742 was obsoleted 5/24/17 and replaced with 1194797 (production effectivity 11/18/16). Service stock of 1128742 can only be used for units built before 11/18/16. The new 1194797 can be used for all.</p> <p>N - 1164808 was replaced with 1194768 (production effectivity 10/13/16). Service stock of 1164808 can only be used for units built before 10/13/16. New 1194768 can be used for all going forward.</p> <p>O - May contain reactivated sieve.</p>				

Quick Reference Chart



When ordering 1130666, you must order 1186230 also.

Common Replacement Filters


Part Number	Description
P0750X030	Bacteria Filter (All Models)
2000489	Cabinet Filter (All Models except Platinum 5)
1107412	Cabinet Filter (Platinum 5 Models and 10 Liter Model)
1143492	Cabinet Filter (Perfecto2/Perfecto2V/PerfectoW)
1131249	Inlet Filter (All Models)

Cabinet, Base, Cylinder and Tabletop


All Models Except Platinum 5




IOH100, IRC6LX02H




IOH250 - Tabletop



**IRC3LX, 5LX, 5LX02, 6LX02,
5LX02T, 6LX02T**



IRC6LX02H



Cabinet, Base, Cylinder and Tabletop

All Models Except Platinum 5

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	2002295	Cabinet Assembly	1	1
		2002001-NLA	NLA - Foam, Side Cabinet	1	1
		2001999-NLA	NLA - Foam, Front Cabinet	1	1
2		2000489	Filter, Cabinet	1	2
3		2002107	Grid, Filter	1	2
4		2001609-NLA	NLA - Filter Grid and Cabinet Filter Assembly	1	2
5		2002540	NLA - Rivet, Plastic	1	1
6		1022500-NLA	NLA - Retainer, Push On	1	1
7		P0875X071-NLA	Strap, Hook and Loop	1	1
8		1070050	Tie Wrap (3-29/32" L)	1	1
9		1008547	Cord, Elastic	1	1
10		2000508	Screw, Phillips Truss Head Self Tap (#8-32 x 5/8")	1	8
11	H	1074004-NLA	NLA - Base Assembly, IRC5LX02T	1	1
11	H	1075661-NLA	NLA - Base Assembly, IRC6LX02	1	1
11	H	2001792-NLA	NLA - Base Assembly, IRC3LX, 5LX, 5LX02	1	1
11	H	1075662-NLA	NLA - Base Assembly, IRC6LX02T	1	1
11	H	1073358-NLA	NLA - Base Assembly, IRC6LX02H	1	1
12	B,C	P0875X500	Caster Assembly (2")	1	4
12	D,J	1076830-NLA	NLA - Caster (2")	1	4
12		1110426	Caster Assembly (2")	1	4
13		2000576	Stand-Off, Twist-Lok (1/2")	1	1
14		2000360	Grommet, Plastic	1	8
15		1056249	Relief, Strain	1	1
16	G	2001761	Cord, Power, Domestic U.S.	1	1
16	F	2000773	Cord, Power, Domestic U.S.	1	1
16		1081182-NLA	NLA - Cord, Power, Canadian	1	1
17		1075469-NLA	NLA - Washer (11/16 x 1/4 x 1/16")	1	1
18		1125090	Label, Oxygen Cylinder	1	1
19	I	IOHM6	Oxygen Cylinder and Regulator, M6 Size Cylinder	1	1
19	I	IOHM9	Oxygen Cylinder and Regulator, C Size Cylinder	1	1
19	I	IOHD	Oxygen Cylinder and Regulator, D Size Cylinder	1	1
19	I	IOHE	Oxygen Cylinder and Regulator, E Size Cylinder	1	1
	E	1079987-NLA	NLA - Cable (Telemetry Set-Up)	1	1
	E	1075248	Label, Cylinder Ring	1	1
	E	1075219	Label, Cylinder	1	1
	E	1013763-NLA	NLA - Label, Power Cord Warning	1	1
20		1101103-NLA	NLA - Top, Table	1	1
21		1102388-NLA	NLA - Leg	1	2
22	K	1105524-NLA	NLA - Brace	1	1
23		1102628-NLA	Kit, Table Top Hardware	1	1
			NLA - Screw, Phillips Pan Head (1/4-20 x 1-1/2")	1	4
			Nut, Wing (1/4-20)	1	4
			NLA - Washer (1/4 x 23/32 x 1/16")	1	4



Cabinet, Base, Cylinder and Tabletop

All Models Except Platinum 5

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	E	1101911-NLA	NLA - Carton, Table	1	1
	E	1100874-NLA	NLA - Instruction Sheet, Table Top	1	1
	E	1056748-NLA	NLA - Suppressor, Surge	1	1
NOTE: A - Includes Items 2 thru 9. B - Attach Casters using Screws. If necessary, order Screw p/n P0851X410 (4 ea) and Keps Nut p/n P0870X400 (4 ea). C - Before 6/15/98 (Screw in Style). D - After 6/14/98 (Stem Style). E - Not Shown. F - To be used on 301/501 Models. G - To be used with LX Models. H - Includes Casters, Foam and Standoff. I - Includes Labels and Carton. J - Item no longer available, refer to Part #1110426. K - Individual Item no longer available, replace with an IOH260 Ready Rack.					

Cabinet and Base

Platinum 5 Models





Cabinet and Base Platinum 5 Models

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	E,A	1108803	Cabinet Cover without Tie Wraps	1	1
1		1115373-NLA	NLA - Cabinet Cover with Tie Wraps	1	1
2	C	1070050	Tie Wrap (3-29/32" L)	1	1
3		1008547	Cord, Elastic	1	1
		1115868	Kit, Cabinet Handle	1	1
4			Handle, Top, Warm Grey	1	1
5			NLA - Screw, Hex Head Tap (1/4-15 x 3/4")	1	2
6		1107412	Filter, Cabinet	1	1
7		1107400	Grid, Filter	1	1
8	E	1107414	Foam, Side Shroud	1	2
9		1107413	Foam, Front Shroud	1	1
10		1172003	Package, Power Cord Strap w/ Hardware	1	1
11			Strap, Hook and Loop	10	1
12			Rivet, Plastic	10	1
			NLA - Retainer, Push On	10	1
13		1110766	Door, Filter Access with Foam Gasket	1	1
14	B	1112258	Base, Domestic	1	1
15		1108816	Foam, Base	1	1
16		1104151	Fitting, Hose Barb with Nut (1/8")	1	1
17		1108211	Cap, End, Dust, Green	1	1
18		1108764	Cord, Power	1	1
19	D	1056249	Relief, Strain	1	1
20		1076830-NLA	NLA - Caster (2")	1	4
20		1110426	Caster Assembly (2")	1	4
21		2000576	Stand-Off, Twist-Lok (1/2")	1	1
22		2000360	Grommet, Plastic	1	4
23	F	2000508	Screw, Phillips Truss Head Self Tap (#8-32 x 5/8")	1	4
		1056748-NLA	NLA - Suppressor, Surge	1	1

NOTE: A - Includes Cover Cabinet and Items 2-12.
 B - Includes Base Cabinet and Items 15-21.
 C - Kit Includes Items 4-5.
 D - Item no longer available, refer to Part Number 1110426.
 E - After 6/28/07, the Hook and Loop Strap are no longer used on the Cabinet Cover Assembly. If Hook and Loop Strap is needed you must order Package Part Number 1172003 for replacement.
 F - Not Shown.

Cabinet and Base

Platinum 10




Cabinet and Base Platinum 10

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,D	1108803	Cabinet Cover without Tie Wraps	1	1
2		1006525	Tie Wrap (3-29/32" L)	1	1
3		1008547	Cord, Elastic	1	1
	C	1115868	Kit, Cabinet Handle	1	1
4			Handle, Top, Warm Grey	1	1
5			NLA - Screw, Hex Head Tap (1/4-15 x 3/4")	1	2
6		1107412	Filter, Cabinet	1	1
7		1107400	Grid, Filter	1	1
8		1107414	Foam, Side Shroud	1	2
9		1107413	Foam, Front Shroud	1	1
10		1172003	Package, Power Cord Strap w/ Hardware	1	1
11			Strap, Hook and Loop	10	1
12			Rivet, Plastic	10	1
			NLA - Retainer, Push On	10	1
13		1110766	Door, Filter Access with Foam Gasket	1	1
14	D	1105499-NLA	NLA - Base, Domestic	1	1
14	D	1112258	Base, Domestic	1	1
15		1108816	Foam, Base	1	1
16		1104151	Fitting, Hose Barb with Nut (1/8")	1	1
17		1108211	Cap, End, Dust, Green	1	1
18		1108764	Cord, Power	1	1
19		1056249	Relief, Strain	1	1
20		1110426	Caster Assembly (2")	1	4
21		2000576	Stand-Off, Twist-Lok (1/2")	1	1
22		2000360	Grommet, Plastic	1	4
23		2000508	Screw, Phillips Truss Head Self Tap (#8-32 x 5/8")	1	4
		1056748-NLA	NLA - Suppressor, Surge	1	1
NOTE: A - Includes Cabinet Cover and Items 2-13. C - After 6/28/07, the Hook and Loop Strap are no longer used on the Cabinet Cover Assembly. If Hook and Loop Strap is needed you must order the Package, Cord Power Strap w/ Hardware Part Number 1172003. D - Includes Base and items 14-21. 1105499 is NLA. Order item #14 part number 1112258 Base Cabinet Assembly as a replacement.					

Cabinet and Base

Perfecto2




Cabinet and Base Perfecto2

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1143470	Top Shroud Assembly, Graphite	1	1
1	A,C,E	1153341	Top Shroud Assembly, Beige	1	1
2		1147920	Label, Text	1	1
3		1137267	Label, Diagnostic	1	1
4		1143492	Filter, Cabinet	1	1
5		1143432-NLA	NLA - Door, Cabinet Access, Graphite	1	1
5	D,E,G	1161831	Door, Cabinet Access, Beige	1	1
6		1143550	Clip, Retaining	1	2
7		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	4
8		1143464	Foam, Shroud - Front	1	1
9		1143465	Foam, Shroud - Side	1	2
10	F	1163095	Base Assembly	1	1
11			Base and Baffle Assembly, Concentrator	1	1
12			Caster Assembly (2")	1	4
13			Cap, End, Dust, Green	1	1
14			Fitting, Hose Barb with Nut (1/8")	1	1
15			NLA - Label, Spec	1	1
16			Cord, Power	1	1
17			Relief, Strain	1	1
11		1168590	Base and Baffle Assembly, Concentrator	1	1
12		1110426	Caster Assembly (2")	1	4
13		1108211	Cap, End, Dust, Green	1	1
14		1104151	Fitting, Hose Barb with Nut (1/8")	1	1
15		1147921-NLA	NLA - Label, Spec	1	1
16		1143493	Cord, Power	1	1
17		1143546	Relief, Strain	1	1
18		2000360	Grommet, Plastic	1	4
	H	2000508	Screw, Phillips Truss Head Self Tap (#8-32 x 5/8")	1	4
NOTE: A - Includes Cabinet and Items 2 thru 9. C - Top Shroud Assembly color was "Graphite" before 1/16/09. D - Use with units "After 10/30/07." E - Top Shroud Assembly color is Beige "After 1/16/2009". F - Includes items 11 thru 17. G - Cabinet Access Door was revised 10/16/09 to allow the Cabinet Filter to be assembled to the outside of the Cabinet Access Door. The Cabinet Access Door is backward compatible in the field for models manufactured before 10/16/09. H - Not shown in illustration. Used to fasten Cabinet to Base.					

Cabinet and Base

Perfecto2V





Cabinet and Base Perfecto2V

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1165488	Top Shroud Assembly	1	1
2		1164905	Label, Text	1	1
3		1137267	Label, Diagnostic	1	1
4		1143492	Filter, Cabinet	1	1
5		1161831	Door, Cabinet Access, Beige	1	1
6		1143550	Clip, Retaining	1	2
7		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	4
8		1143464	Foam, Shroud - Front	1	1
9		1143465	Foam, Shroud - Side	1	2
10	B	1163095	Base Assembly	1	1
11			Base and Baffle Assembly, Concentrator	1	1
12			Caster Assembly (2")	1	4
13			Cap, End, Dust, Green	1	1
14			Fitting, Hose Barb with Nut (1/8")	1	1
15			NLA - Label, Spec	1	1
16			Cord, Power	1	1
17			Relief, Strain	1	1
11		1168950	Base and Baffle Assembly, Concentrator	1	1
12		1110426	Caster Assembly (2")	1	4
13		1108211	Cap, End, Dust, Green	1	1
14		1104151	Fitting, Hose Barb with Nut (1/8")	1	1
15		1147921-NLA	NLA - Label, Spec	1	1
16		1143493	Cord, Power	1	1
17		1143546	Relief, Strain	1	1
18		2000360	Grommet, Plastic	1	4
NOTE: A - Includes Cabinet and Items 2 thru 9. B - Includes items 11 thru 17.					

Cabinet and Base

PerfectoW




Cabinet and Base PerfectoW

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	F,A	1161717-NLA	NLA - Top Shroud Assembly, Beige	1	1
2		1161721-NLA	NLA - Label, Text, English - PerfectoW	1	1
3		1137267	Label, Diagnostic	1	1
4		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	4
5		1153695	Foam, Inlet - Ceiling	1	1
6		1143492	Filter, Cabinet	1	1
7	F,G	1161831	Door, Cabinet Access, Beige	1	1
8		1153700	Foam, Inlet - Wall	1	1
9		1143465	Foam, Shroud - Side	1	2
10		1143464	Foam, Shroud - Front	1	1
11	B E C	1143550	Clip, Retaining	1	2
12		1156532	Inlet Air Scoop Assembly	1	1
		1161718-NLA	NLA - Instruction Sheet, Inlet Air Scoop Assembly	1	1
13		1143429-NLA	NLA - Base, Concentrator	1	1
14		1153706	Foam, Damping - Base	1	2
15		1154501	Foam, Resonator	1	1
16		1153705	Foam, Damping Outlet Base	1	1
17		2000360	Grommet, Plastic	1	4
18		1108211	Cap, End, Dust, Green	1	1
19		1147921-NLA	NLA - Label, Spec	1	1
20		1104151	Fitting, Hose Barb with Nut (1/8")	1	1
21		1110426	Caster Assembly (2")	1	4
22		1153704	Damper, Sound Base	1	1
23		1143444-NLA	NLA - Tech Only - Baffle, Exhaust	1	1
24	1143546	Relief, Strain	1	1	
25	1143493	Cord, Power	1	1	
NOTE: A - Includes Cabinet and Items 2-11. B - Includes Inlet Foam and Scoop. When ordering item 12, order instruction sheet part number 1161718 also. C - Includes Base Only. E - Not Shown. F - Changed Top Shroud Assembly Color back to Beige. After 1/16/2009. G - Cabinet Access Door was revised 10/16/09 to allow the Cabinet Filter to be assembled to the outside of the Cabinet Access Door. The Cabinet Access Door is backward compatible in the field for models manufactured before 10/16/09.					


Compressor and Mounting Hardware

All Models Except Platinum 5


3LX (After 7/10/97)
5LX (After 4/21/97)
6LX (All Dates)



3LX (Before 7/11/97)



5LX (Before 4/22/97)



SRE_005A

Compressor and Mounting Hardware

All Models Except Platinum 5

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1075071-NLA	NLA - Compressor Assembly - 6LX, 6LX02	1	1
1		1075071E-NLA	NLA - WARR RPLMT - Compressor Assembly (6LX, 6LX02)	1	1
	A,F	2000894-NLA	NLA - Kit, Compressor Rebuild	1	1
			Kit, Rebuild	1	1
			NLA - Instruction Sheet	1	1
	J	1157138	Kit, Compressor Complete Assembly - Thomas	1	1
			Compressor with Mounting Hardware - Thomas	1	1
			NLA - Tech Only - Capacitor, 15 MFD	1	1
			Tie Wrap (8" L)	2	1
		1186230	Kit, Compressor Restraint Bracket w/ Hardware	1	1
			Bracket, Restraint, Compressor	1	1
			NLA - Bumper, Bracket	1	1
			Bumper, Restraint Bracket, Compressor	1	1
			NLA - Screw, Phillips Hex Washer Head w/ Patch (1/4-20 x 3/4")	1	1
			Retainer, Push (1/4")	1	1
		1106690E-NLA	NLA - WARR RPLMT - Compressor Assembly - Thomas (5LX Plat.)	1	1
	A,G,I	1109235-NLA	Kit, Rebuild for Thomas Compressors	1	1
			Kit, Rebuild	1	1
			NLA - Instruction Sheet	1	1
	A,G	1154510	Kit, Compressor Rebuild w/ Trigon Seals & Instruction Sheet - IRC5LX/IRC5LX02 (Thomas)	1	1
			Compressor Rebuild Kit w/ Trigon Seals - IRC5LX/IRC5LX02 (Thomas)	1	1
			Instruction Sheet, Thomas Rebuild 2660-2650	1	1
	A,H	2000904-NLA	NLA - Kit, Compressor Rebuild	1	1
			NLA - Kit, Rebuild	1	1
			NLA - Instruction Sheet	1	1
2		1086187-NLA	NLA - Hose, Barbed, Patch Lock (1/2 x 1/4" NPT)	1	1
3	B	2001374-NLA	NLA - Fitting, Resonator	1	1
3	C	1108819	Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	1
4		P0876X026	Clamp, Hose	1	1
5	B	1095537-NLA	NLA - Kit, Resonator Retro-Fit	1	1
			Resonator	1	1
			NLA - Screw, Phillips Pan Head (#10-16 x 1/2")	1	2
			NLA - Tubing, Silicon (1/2 x 5")	1	1
			Tube, Braided (10 x 1/2")	1	1
			Spring, Extension (7/16 x 9-3/4" L)	1	1
			Screw, Socket Button Head (1/4-20 x 3/8")	1	1
			NLA - Bracket, Mounting	1	1
			Clip, Flange	1	1
			Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	1

Compressor and Mounting Hardware

All Models Except Platinum 5


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
			Tie Wrap (8" L)	1	1
			NLA - Instruction Sheet	1	1
5	C	1085873	Resonator	1	1
6		2001682	Tie Wrap (8" L)	1	3
7		2002264	Tubing, Reinforced Braided (1/2 x 26")	1	1
8		1088316	Spring, Extension (7/16 x 9-3/4" L)	1	1
9		2000888-NLA	NLA - Mounting Bracket Assembly	1	2
10		2000876-NLA	NLA - Washer, Cup	1	4
11		1105418-NLA	NLA - Screw, Phillips Pan Head (1/4-20 x 1-1/2")	1	4
12		1085696	NLA - Screw, Phillips Truss Head (1/4-20 x 1/2")	1	4
13		P0875X510	Mount, Diamond	1	4
14		2002906-NLA	NLA - Spacer (1/4 x 1/2 x 7/32")	1	4
15		P0875X560-NLA	NLA - Boot, Terminal	1	1
16	E	2000770-NLA	NLA - Capacitor - MFD3LX	1	1
16	D	2003019	Kit, Capacitor 15 MFD	1	1
18			NLA - Tech Only - Capacitor, 15 MFD	1	1
19			Clamp, Cap	1	1
20			Screw, Phillips Pan Head (8-32 x 1/2")	1	1
21			Washer, Flat (#8)	1	1
22			NLA - Nut, Keps (#6-32)	1	1
23			NLA - Boot, Terminal	1	1
24			Tie Wrap (8" L)	1	2
25			Tie Wrap (3-29/32" L)	1	1
26			Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	1
17	B	1085697-NLA	NLA - Relief Valve w/ Patch Lock, 35 PSIG	1	1

NOTE: A - Not Shown.
 B - Before 7/00.
 C - After 6/00.
 D - Use with IRC5LX/ IRC6LX After 1/97 and IRC3LX After 6/97.
 E - Use with IRC5LX/ IRC6LX Before 2/97 and IRC3LX Before 7/97.
 F - Rebuild Kit is for Thomas Compressor Model #2618 and #2638. Compressor Model numbers are required when ordering the Rebuild Kit.
 G - Rebuild Kit is for Thomas Compressor Model #2650 and #2660. Compressor Model numbers are required when ordering the Rebuild Kit.
 H - Rebuild Kit is for Thomas Compressor Model #2619 and #2639. Compressor Model numbers are required when ordering the Rebuild Kit.
 I - Part #1109235 is No Longer Available. Order part #1154510 as a replacement kit.
 J - Includes Compressor, Capacitor, Tie Wraps, and 1186230.

Compressor and Mounting Hardware

Platinum 5 Models

NOTE: Compressors are interchangeable. However, Capacitor and Bracket Restraint must also be replaced at the same time.



SRE_013A

Compressor and Mounting Hardware

Platinum 5 Models

NOTE: Compressors are interchangeable. However, Capacitor and Bracket Restraint must also be replaced at the same time.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1157138	Kit, Compressor Complete Assembly - Thomas	1	1
			Compressor with Mounting Hardware - Thomas	1	1
			NLA - Tech Only - Capacitor, 15 MFD	1	1
			Tie Wrap (8" L)	2	1
		1186230	Kit, Compressor Restraint Bracket w/ Hardware	1	1
			Bracket, Restraint, Compressor	1	1
			NLA - Bumper, Bracket	1	1
			Bumper, Restraint Bracket, Compressor	1	1
			NLA - Screw, Phillips Hex Washer Head w/ Patch (1/4-20 x 3/4")	1	1
			Retainer, Push (1/4")	1	1
1	B	1106690E-NLA	NLA - WARR RPLMT - Compressor Assembly - Thomas (5LX Plat.)	1	1
2		1108819	Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	2
3		1105504-NLA	NLA - Spring, Compressor Mount	1	4
3		1150581	Spring, Compressor Mount	1	4
4		1106446-NLA	NLA - Pin, Spring Retainer	1	4
4		1150580	Pin, Spring Retainer	1	4
5		1107411	Glide, Pin	1	4
6		1186230	Kit, Compressor Restraint Bracket w/ Hardware	1	1
7			Bracket, Restraint, Compressor	1	1
8			NLA - Bumper, Bracket	1	1
9			Bumper, Restraint Bracket, Compressor	1	1
10			NLA - Screw, Phillips Hex Washer Head w/ Patch (1/4-20 x 3/4")	1	1
11			Retainer, Push (1/4")	1	1
	C,D,B	1154510	Kit, Compressor Rebuild w/ Trigon Seals & Instruction Sheet - IRC5LX/IRC5LXO2 (Thomas)	1	1
			Compressor Rebuild Kit w/ Trigon Seals - IRC5LX/IRC5LXO2 (Thomas)	1	1
			Instruction Sheet, Thomas Rebuild 2660-2650	1	1
	D,E	1068698-NLA	NLA - Tech Only - Capacitor, 15 MFD	1	1
	D	2001682	Tie Wrap (8" L)	1	2
12	F	1125371-NLA	NLA - Compressor with Mounting Hardware - Gast	1	1
12	G	1125371E-NLA	NLA - Compressor Assembly, Warranty Replacement - Gast (5LX Plat.)	1	1
	D,H	1124028	Kit, Rebuild for Gast Compressors	1	1
			Kit, Rebuild	1	1
			Instruction Sheet - Compressor Rebuild	1	1
	D,I,J	1131250	Kit, Capacitor (10 MFD)	1	1
			NLA - Tech Only - Capacitor, 10 MFD	1	1
			Tie Wrap (17-1/2" L)	1	1

Compressor and Mounting Hardware

Platinum 5 Models


NOTE: Compressors are interchangeable. However, Capacitor and Bracket Restraint must also be replaced at the same time.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
13	J	1121244	Tie Wrap (8" L)	1	2
			Instruction Sheet, Capacitor Replacement	1	1
			Kit, Compressor Restraint Bracket (Gast)	1	1
NOTE: A - Kit includes items 2 - 6 and 1186230. Must use with Thomas Compressors. B - When replacing or servicing a Compressor, the Capacitor (1068698) and two Tie Wraps (2001682) must also be replaced. C - Kit includes Retainer Plates and Trigon Seals to enable rebuilding of both Trigon and pre-Trigon configurations. D - Not shown. E - Must use with Thomas Compressors. F - Item is No Longer Available, refer to item 1 (1157138) as its replacement. G - When replacing or servicing a Compressor, the Capacitor (1131250) must also be replaced. H - Includes Instruction Sheet, (2) Head O-Rings, (4) Leaf Valves, (4) Self Tapping Screws, (4) Retainer Screws, and (2) Piston Cups. I - Includes Capacitor, Tie Wraps, and Installation Instructions. J - Must use with Gast Compressors.					

Compressor and Mounting Hardware

Platinum 10

When ordering 1130666 or 1130666E, you must order Kit part number 1186230.



Compressor and Mounting Hardware

Platinum 10

When ordering 1130666 or 1130666E, you must order Kit part number 1186230.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1130666	Compressor Assembly - 10LX/ 10LXO2	1	1
1	E	1130666E	WARR RPLMT - Compressor Assembly - 10LX/ 10LXO2	1	1
	B,C	1153687	Kit, Compressor Rebuild w/ Instructions (2660CE50)	1	1
			Compressor Rebuild Assembly - Thomas	1	1
			Instruction Sheet, Thomas Rebuild 2660-2650	1	1
	C	1128740	Capacitor, 30 MFD	1	1
	C	1089455	Tie Wrap (15-1/4" L)	1	2
	C,G	1188159	Kit, Compressor Rebuild w/ Instructions (ZW400D2)	1	1
			Manual, Sevice (ZW400D2)	1	1
			Cylinder	1	2
			Seal, Piston Cup	1	2
			O-Ring, Head Seal	1	2
			Screw, Valve Plate	1	4
			Washer, Square	1	4
			Washer, Semi-Circle Reinforced Metal	1	2
			Washer, Plastic Semi-Circle Reinforced	1	2
			Washer, Metal Square Reinforced	1	2
			Flapper, Intake Valve (Narrow)	1	2
			Flapper, Exhaust Valve (Wide)	1	2
			O-Ring, Valve Plate	1	2
			Screw, Head	1	8
	C	1128740	Capacitor, 30 MFD	1	1
	C	1089455	Tie Wrap (15-1/4" L)	1	2
2		1128741	Compressor - 10LX/10LXO2	1	1
2		1128741E-NLA	NLA - WARR RPLMT - Compressor - 10LX/10LXO2	1	1
3		1108819	Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	2
4		1105504-NLA	NLA - Spring, Compressor Mount	1	4
4	H	1150581	Spring, Compressor Mount	1	4
5		1106446-NLA	NLA - Pin, Spring Retainer	1	4
5	H	1150580	Pin, Spring Retainer	1	4
6		1107411	Glide, Pin	1	4
7	D,F	1186230	Kit, Compressor Restraint Bracket w/ Hardware	1	1
8			Bracket, Restraint, Compressor	1	1
9			NLA - Bumper, Bracket	1	1
10			Bumper, Restraint Bracket, Compressor	1	1
11			NLA - Screw, Phillips Hex Washer Head w/ Patch (1/4-20 x 3/4")	1	1
12			Retainer, Push (1/4")	1	1

Compressor and Mounting Hardware


Platinum 10

When ordering 1130666 or 1130666E, you must order Kit part number 1186230.

NOTE: A - Includes item 2 thru item 7.
B - When replacing or servicing a Compressor, the Capacitor (1128740) and two Tie Wraps (1089455) must also be replaced.
Used for Thomas Compressors only.
C - Not Shown
D - Kit includes items 8 - 12.
E - When ordering 1130666E, you must order Kit part number 1186230.
F - Included with Item 1 (1130666) only.
G - When replacing or servicing a Compressor, the Capacitor (1128740) and two Tie Wraps (1089455) must also be replaced.
Used for GSE Compressors (ZW400D2) only.
H - It is recommended to order item 4 and 5 together. The newer Spring and Pin is now narrower.

Compressor and Mounting Hardware

Perfecto2/PerfectoW



Compressor and Mounting Hardware


Perfecto2/PerfectoW

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,D	1162038	Kit, Compressor with Hardware - PW/PO2W	1	1
			NLA - Compressor - P2/ PO2/PW/P02W	1	1
2			Screw, Hex Head Flange (1/4-20 x 7/8")	1	4
3			Bushing, Mounting	1	4
4			Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	1
5			Elbow, Male Flare	1	1
6			Spring, Compression (63/64 x 1-1/4 x 1/16")	1	4
7			NLA - Tech Only - Capacitor, 10 MFD	1	1
8	A	1162038E	Tie Wrap (17-1/2" L)	1	1
1			WARR RPLMT - Kit, Compressor w/ Hardware, PW/PO2W	1	1
			NLA - WARR RPLMT - Compressor, P2/PO2/PW/PO2W	1	1
2			Screw, Hex Head Flange (1/4-20 x 7/8")	1	4
3			Bushing, Mounting	1	4
4			Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	1
5			Elbow, Male Flare	1	1
6			Spring, Compression (63/64 x 1-1/4 x 1/16")	1	4
7	E	1169408-NLA	NLA - Tech Only - Capacitor, 10 MFD	1	1
8			Tie Wrap (17-1/2" L)	1	1
2			Screw, Hex Head Flange (1/4-20 x 7/8")	1	4
3			Bushing, Mounting	1	4
4			Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	1
5			NLA - Fitting, Elbow Swivel (1/4 NPT x 3/8 Tube)	1	1
5			Elbow, Male Flare	1	1
6			Spring, Compression (63/64 x 1-1/4 x 1/16")	1	4
7	B	1131250	Kit, Capacitor (10 MFD)	1	1
			NLA - Tech Only - Capacitor, 10 MFD	1	1
			Tie Wrap (17-1/2" L)	1	1
			Tie Wrap (8" L)	1	2
			Instruction Sheet, Capacitor Replacement	1	1
8			Tie Wrap (17-1/2" L)	1	1
			Kit, Compressor Rebuild w/ Instruction Sheet	1	1
			Compressor Rebuild - Thomas	1	1
	C	1153688	Instruction Sheet, Thomas Rebuild 2660-2650	1	1

NOTE: A - Kit includes the Compressor and items 2 thru 8. Compressor is not available separately.
 B - Includes 10MFD Capacitor, Tie Wraps and Installation Instructions.
 C - When replacing or servicing a compressor, Part Number 1131250 must also be replaced.
 D - When ordering a replacement compressor, you should order 1143557 Heat Exchanger when replacing a unit manufactured between 5/17/11 and 7/20/11.
 E - No longer available, use 1143489 as its replacement.

Compressor and Mounting Hardware

Perfecto2V




Compressor and Mounting Hardware Perfecto2V

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,D	1165097	Kit, Compressor with Hardware - PO2V	1	1
			Compressor - PO2V	1	1
2			Screw, Hex Head Flange (1/4-20 x 7/8")	1	4
3			Bushing, Mounting	1	4
4			Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	1
5			Elbow, Male Flare	1	1
6			Spring, Compression (63/64 x 1-1/4 x 1/16")	1	4
7			Capacitor, 22MFD	1	1
8	E	1169408-NLA	Tie Wrap (17-1/2" L)	1	1
2			Screw, Hex Head Flange (1/4-20 x 7/8")	1	4
3			Bushing, Mounting	1	4
4			Fitting, Elbow (1/2" Beaded Hose x 1/4" NPT)	1	1
5			NLA - Fitting, Elbow Swivel (1/4 NPT x 3/8 Tube)	1	1
5			Elbow, Male Flare	1	1
6			Spring, Compression (63/64 x 1-1/4 x 1/16")	1	4
9	B	1191077	Kit, Capacitor Replacement, 22MFD	1	1
			Capacitor, 22MFD	1	1
			Tie Wrap (17-1/2" L)	1	1
			Tie Wrap (8" L)	1	2
			Instruction Sheet, Capacitor Replacement	1	1
			0	1	1
			Kit, Cylinder Rebuild, Complete, PO2V	1	1
			Kit, Cylinder w/ Hardware, PO2V	1	1
	F	1180335	Package, O-Ring, PO2V	100	1


NOTE: A - Kit includes the Compressor and items 2-8
B - Kit includes items 7 and 8
C - Not Shown. Kit includes (2) ea. Cylinders, (2) ea. Seal Cups, (4) ea. Reed Valves, (4) ea. O-Rings, (12) ea. Head Mounting Screws, (4) ea. Valve Flapper Screws, (4) Washers and Instruction Sheet
D - When ordering a replacement compressor, you should order 1143557 Heat Exchanger, when replacing a unit manufactured between 5/17/11 and 7/20/11.
E - No longer available. Use 1143489 as its replacement.
F - Not shown in illustration.

Control Panel and P.C. Board

All Models Except Platinum 5 (Before 12/15/00)




1128379



Used with Flowmeters that
have an Indented Phillips
Head Screw on top.

1128380



Used with Flowmeters that
have a Slotted Head Screw
on top.

Control Panel and P.C. Board

All Models Except Platinum 5 (Before 12/15/00)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1101585-NLA	NLA - Control Panel Assembly - IRC3LX	1	1
1		1101584	Control Panel Assembly - IRC5LX	1	1
1		1101583	Control Panel Assembly - IRC5LX02	1	1
1		1073754-NLA	NLA - Control Panel Assembly - IRC5LX02T	1	1
1		1101582	Control Panel Assembly - IRC6LX02/ IRC6LX02H	1	1
1		1073752	Control Panel Assembly - IRC6LX02T	1	1
2		2001769-NLA	NLA - Panel, Control	1	1
3		40120X001	Screw, Phillips Pan Head (#8-32 x 3/8")	1	2
4		1080809-NLA	Label, Control Panel, English - IRC3LX	1	1
4		2002356-NLA	NLA - Label, Control Panel, English - IRC5LX	1	1
4		2002355-NLA	NLA - Label, Control Panel, English - IRC5LX02, 5LX02T	1	1
4		1074670	Label, Control Panel, English - IRC6LX02, 6LX02T, 6LX02H	1	1
4		1079329	Label, Control Panel, French Canadian - IRC5LX	1	1
4		1079330	Label, Control Panel, French Canadian - IRC5LX02	1	1
4		1079331	Label, Control Panel, French Canadian - IRC6LX02	1	1
4		1079324	Label, Control Panel, Spanish - IRC5LX	1	1
4		1079325	Label, Control Panel, Spanish - IRC5LX02	1	1
4		1079326	Label, Control Panel, Spanish - IRC6LX02	1	1
5		1180107	Connector, Barbed, Brass (1/4")	1	1
6		2001086-NLA	NLA - Knob, Flowmeter	1	1
6	I	1128379	Package, Knob, Flowmeter, Silver w/ Black Swish (Pkg. of 10)	10	1
6	J	1128380	Package, Knob, Flowmeter, Black w/ Silver Swish	10	1
7		2001989-NLA	NLA - Flowmeter with Palnuts - 3 Liter	1	1
7	H	2001773-NLA	NLA - Flowmeter with Palnuts - 5 Liter	1	1
7	G	1110561	Flowmeter without Palnuts - 5 Liter	1	1
7		1074672-NLA	NLA - Flowmeter with Palnuts - 6 Liter	1	1
8		2000799	Adapter with Tubing for Humidifier Bottle	1	1
9		2000503	Adapter, Humidifier Bottle	1	1
10		1098707	Hour Meter Assembly with Clip	1	1
11		2001134	Switch, On/Off	1	1
12		1056554	Circuit Breaker (8 AMP)	1	1
13		2000464-NLA	NLA - Support, Control Panel	1	2
14		2002674	Palnut (3/8-27)	1	1
15		2000507	Palnut (9/16-18)	1	2
16		1061773-NLA	NLA - Remote LED Assembly, Sens O2 - IRC5LX02	1	1
17	B	2001696-NLA	NLA - Transformer Assembly	1	1
17	C	1085200	Transformer Assembly	1	1
18		1028794-NLA	NLA - Screw, Phillips Truss Head Tap (#6 x 1/4")	1	6
19		2000436	Check Valve, Barb (5/16")	1	2
20		2001682	Tie Wrap (8" L)	1	2

Control Panel and P.C. Board

All Models Except Platinum 5 (Before 12/15/00)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
21		P0750X030	Filter, HEPA Outlet	1	1
22		2000897	Tubing, Tygon (1/4 x 24")	1	1
23	D	2001834-NLA	NLA - Harness, Wire IRC5LX	1	1
23	D	2001762-NLA	NLA - Harness, Wire-Sensor - IRC5LX02	1	1
23	E	2003085-NLA	NLA - Harness, Wire, (All LX Models)	1	1
24		1036162	Stand-Off, P.C. Board	1	6
25	F	See Chart	# 25 P.C. Boards	1	1
	A	2001321-NLA	NLA - T-Fitting	1	1

NOTE: A - Not Shown.
 B - Before 11/97.
 C - After 10/97.
 D - Before 5/97.
 E - After 4/97.
 F - See P.C. Board Chart.
 G - Palnuts are no longer included, must be ordered separately. Requires (2) Palnuts - #2000507.
 H - Item no longer available, refer to Part #1110561.
 I - Knob is silver with black "swish". Used on Flowmeters with a Phillips Head Indented Screw on top.
 J - Knob is black with silver "swish". Used on Flowmeters with a Slotted Head Screw on top.


#25 P.C. Boards

Model	P.C. Board	Reconditioned P.C. Board
IRC3LX	1078183-NLA	1078183E
IRC5LX	1078183-NLA	1078183E
IRC5LX02 Before 7/98	1078181-NLA	1078181E-NLA
IRC5LX02 After 6/98	1113693	1113693E-NLA
IRC5LX02T	1079444-NLA	N/A
IRC6LX02 Before 7/98	1078181-NLA	1078181E-NLA
IRC6LX02 After 6/98	1113693	1113693E-NLA
IRC6LX02T	1079444-NLA	N/A

All P.C. Boards include an Instruction Sheet.

Control Panel and P.C. Board

All Models Except Platinum 5 (After 12/14/00)




1128379



Used with Flowmeters that
have an Indented Phillips
Head Screw on top.

1128380



Used with Flowmeters that
have a Slotted Head Screw
on top.

Control Panel and P.C. Board

All Models Except Platinum 5 (After 12/14/00)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1101585-NLA	NLA - Control Panel Assembly - IRC3LX	1	1
1		1101584	Control Panel Assembly - IRC5LX	1	1
1		1101583	Control Panel Assembly - IRC5LX02	1	1
1		1101582	Control Panel Assembly - IRC6LX02/ IRC6LX02H	1	1
2		1081276	Control Panel Only	1	1
3		2003065-NLA	NLA - Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
4		1180107	Connector, Barbed, Brass (1/4")	1	1
5		2001086-NLA	NLA - Knob, Flowmeter	1	1
5	D	1128379	Package, Knob, Flowmeter, Silver w/ Black Swish (Pkg. of 10)	10	1
5	E	1128380	Package, Knob, Flowmeter, Black w/ Silver Swish	10	1
6		2000799	Adapter with Tubing for Humidifier Bottle	1	1
7		2000503	Adapter, Humidifier Bottle	1	1
8		2001989-NLA	NLA - Flowmeter with Palnuts - 3 Liter	1	1
8	C	2001773-NLA	NLA - Flowmeter with Palnuts - 5 Liter	1	1
8	B	1110561	Flowmeter without Palnuts - 5 Liter	1	1
8		1074672-NLA	NLA - Flowmeter with Palnuts - 6 Liter	1	1
9		1098707	Hour Meter Assembly with Clip	1	1
10		2001134	Switch, On/Off	1	1
11		1101587-NLA	NLA - Label, Text Control Panel - IRC3LX	1	1
11		1105530-NLA	NLA - Label, Text Control Panel - IRC5LX	1	1
11		1105532-NLA	NLA - Label, Text Control Panel - IRC5LX02	1	1
11		1101586-NLA	NLA - Label, Text Control Panel - IRC6LX02	1	1
12		1101590-NLA	NLA - Label, Logo Control Panel - IRC3LX	1	1
12		1105531-NLA	NLA - Label, Logo Control Panel - IRC5LX, 5LX02	1	1
12		1101588-NLA	NLA - Label, Logo Control Panel - IRC6LX02	1	1
13		1101591	Support, Control Panel	1	2
14		1056554	Circuit Breaker (8 AMP)	1	1
15		2000507	Palnut (9/16-18)	1	2
16		2002674	Palnut (3/8-27)	1	1
17		1085200	Transformer Assembly	1	1
18		1102696	Screw, Phillips Pan Head SEMS (#8-32 x 3/8")	1	2
19		1101620-NLA	NLA - Harness, Wire (24V)	1	1
20		P0750X030	Filter, HEPA Outlet	1	1
21		2000897	Tubing, Tygon (1/4 x 24")	1	1
22		2001682	Tie Wrap (8" L)	1	2
23	A	See Chart	#23 P.C. Boards	1	1

NOTE: A - See P.C. Board Chart.
B - Palnuts not included. See item 15.
C - Item no longer available, refer to Part #1110561.
D - Knob is silver with black "swish". Used on Flowmeters with an Indented Phillips Head Screw on top.
E - Knob is black with silver "swish". Used on Flowmeters with a Slotted Head Screw on top.



Control Panel and P.C. Board

All Models Except Platinum 5 (After 12/14/00)

#23 P.C. Boards

Model	P.C. Board	Reconditioned P.C. Board
IRC3LX	1078183-NLA	1078183E
IRC5LX	1078183-NLA	1078183E
IRC5LX02	1113693	1113693E-NLA
IRC5LX02T	1079444-NLA	N/A
IRC6LX02	1113693	1113693E-NLA
IRC6LX02T	1079444-NLA	N/A
IRC6LX02H	1113693	1113693E-NLA
All P.C. Boards include an Instruction Sheet.		

Control Panel and P.C. Board (Platinum 5 Models) Before 12/5/02




1128379



Used with Flowmeters that have an Indented Phillips Head Screw on top.

1128380



Used with Flowmeters that have a Slotted Head Screw on top.

Control Panel and P.C. Board

(Platinum 5 Models)

Before 12/5/02

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	J	1110784-NLA	NLA - Control Panel Assembly - IRC5LX	1	1
1		1115189-NLA	NLA - Control Panel Assembly - IRC5LX	1	1
1	I	1106693-NLA	NLA - Control Panel Assembly - IRC5LX02	1	1
1	I	1115190-NLA	NLA - Control Panel Assembly - IRC5LX02	1	1
2		2000507	Palnut (9/16-18)	1	3
3		1056554	Circuit Breaker (8 AMP)	1	1
4		1098707	Hour Meter Assembly with Clip	1	1
5		2001134	Switch, On/Off	1	1
6		1105530-NLA	NLA - Label, Text Control Panel - IRC5LX	1	1
6		1105532-NLA	NLA - Label, Text Control Panel - IRC5LX02	1	1
7	F	2001773-NLA	NLA - Flowmeter with Palnuts - 5 Liter	1	1
7		1110561	Flowmeter without Palnuts - 5 Liter	1	1
8	G	1128379	Package, Knob, Flowmeter, Silver w/ Black Swish (Pkg. of 10)	10	1
8	H	1128380	Package, Knob, Flowmeter, Black w/ Silver Swish	10	1
9		1180107	Connector, Barbed, Brass (1/4")	1	1
10	D,N	See Note	See Footnote	1	1
10	E	1123239-NLA	NLA - Decal, Logo, Control Panel	1	1
11		2003065-NLA	NLA - Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
11		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
12		1085200	Transformer Assembly	1	1
13		1102696	Screw, Phillips Pan Head SEMS (#8-32 x 3/8")	1	2
14		2001902-NLA	NLA - Screw, Phillips Pan Head (#4-40 x 1/4")	1	2
15		1076275	Tubing, Tygon (1/8 x 18")	1	1
16		1110481-NLA	NLA - Harness, Wire, 9 Pin	1	1
17	B,C	1078183-NLA	NLA - P.C. Board Replacement LX	1	1
17	B,C	1113693	P.C. Board, LX02	1	1
18	A,M	1110779-NLA	NLA - Tubing Assembly, Check Valve Discharge	1	1
19	M	1108815-NLA	NLA - Tubing, Silicon (1/4 x 16-1/4")	1	1
20		1108809	Tubing, PVC (1/4 x 12-3/4")	1	1
21	K,L	2001682	Tie Wrap (8" L)	1	2
22	K,L	1110473	Tubing, PVC (1/4 x 1-3/4")	1	2
23	K,L	2000436	Check Valve, Barb (5/16")	1	1
24	K,L	P0750X030	Filter, HEPA Outlet	1	1
25	K,L	1143144	Tubing, PVC (1/4 x 5")	1	1
26		1081276	Control Panel Only	1	1



Control Panel and P.C. Board

(Platinum 5 Models)

Before 12/5/02

NOTE: A - Includes (2) Tubing, Check Valve, Elbow and Tie Wrap.
B - Includes Item 15 and Instruction Sheet.
C - Before 11/19/02.
D - Before 4/14/04.
E - After 4/13/04.
F - Item no longer available, refer to Part #1110561. Palnuts must be ordered separately, refer to Part #2000507 in a quantity of 2.
G - Knob is silver with black "swish". Used on Flowmeters with an Indented Phillips Head Screw on top.
H - Knob is black with silver "swish". Used on Flowmeters with a Slotted Head Screw on top.
I - Item no longer available.
J - Item no longer available, refer to Part #1115189.
K - "After 7/19/06" the location of the HEPA Filter has changed. The HEPA Filter is now located inside the cabinet.
L - For use with Control Panel Assembly "After 7/19/06."
M - For use with Control Panel Assembly "Before 7/20/06."
N - No longer available as a service part in the field.

Control Panel and P.C. Board (Platinum 5 Models) After 12/5/02




1128379



Used with Flowmeters that have an Indented Phillips Head Screw on top.

1128380



Used with Flowmeters that have a Slotted Head Screw on top.

Control Panel and P.C. Board

(Platinum 5 Models)

After 12/5/02

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C	1115189-NLA	NLA - Control Panel Assembly - IRC5LX	1	1
1	D	1115190-NLA	NLA - Control Panel Assembly - IRC5LX02	1	1
2		2000507	Palnut (9/16-18)	1	3
3		1056554	Circuit Breaker (8 AMP)	1	1
4		1098707	Hour Meter Assembly with Clip	1	1
5		2001134	Switch, On/Off	1	1
6		1108820-NLA	NLA - Label, Text Control Panel - IRC5LX	1	1
6		1108755-NLA	NLA - Label, Text Control Panel - IRC5LX02	1	1
7	B	2001773-NLA	NLA - Flowmeter with Palnuts - 5 Liter	1	1
7		1110561	Flowmeter without Palnuts - 5 Liter	1	1
8	H	1128379	Package, Knob, Flowmeter, Silver w/ Black Swish (Pkg. of 10)	10	1
8	G	1128380	Package, Knob, Flowmeter, Black w/ Silver Swish	10	1
9		1180107	Connector, Barbed, Brass (1/4")	1	1
10	F	1123239-NLA	NLA - Decal, Logo, Control Panel	1	1
11		2003065-NLA	NLA - Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
11		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
12		1085200	Transformer Assembly	1	1
13		1102696	Screw, Phillips Pan Head SEMS (#8-32 x 3/8")	1	2
14		2001902-NLA	NLA - Screw, Phillips Pan Head (#4-40 x 1/4")	1	2
15		1076275	Tubing, Tygon (1/8 x 18")	1	1
16		1108737	Harness, Wire, 10 Pin	1	1
17	C	1108720-NLA	NLA - LXN PCB Assembly	1	1
17	C	1108720E-NLA	NLA - WARR RPLMT - LXN PCB Assembly	1	1
17	D	1108731	ISBN PCB Assembly	1	1
17	D	1108731E-NLA	NLA - WARR RPLMT - ISBN PCB Assembly	1	1
18	A,K	1110779-NLA	NLA - Tubing Assembly, Check Valve Discharge	1	1
19	K	1108815-NLA	NLA - Tubing, Silicon (1/4 x 16-1/4")	1	1
20		1108809	Tubing, PVC (1/4 x 12-3/4")	1	1
21	I,J	2001682	Tie Wrap (8" L)	1	2
22	I,J	1110473	Tubing, PVC (1/4 x 1-3/4")	1	2
23	I,J	2000436	Check Valve, Barb (5/16")	1	1
24	I,J	P0750X030	Filter, HEPA Outlet	1	1
25	I,J	1143144	Tubing, PVC (1/4 x 5")	1	1
26		1081276	Control Panel Only	1	1



Control Panel and P.C. Board

(Platinum 5 Models)

After 12/5/02

NOTE: A - Includes (2) Tubing, Check Valve, Elbow and Tie Wrap.
B - Item no longer available, refer to Part #1110561.
C - To be used on model IRC5LX. Part #1108720 is NLA. See part #1108720E as a replacement.
D - To be used on model IRC5LX02.
F - After 4/13/04.
G - Knob is black with silver "swish". Used on Flowmeters with a Slotted Head Screw on top.
H - Knob is silver with black "swish". Used on Flowmeters with an Indented Phillips Head Screw on top.
I - "After 7/19/06" the location of the HEPA Filter has changed. The HEPA Filter is now located inside the cabinet.
J - For use with Control Panel Assembly "After 7/19/06."
K - For use with Control Panel Assembly "Before 7/20/06."

Control Panel and P.C. Board Platinum 10



SRE_021A

Control Panel and P.C. Board

Platinum 10

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1132056-NLA	NLA - Control Panel Assembly - 10LX	1	1
1		1130682	Control Panel Assembly - 10LXO2	1	1
2		2000507	Palnut (9/16-18)	1	3
3		1056554	Circuit Breaker (8 AMP)	1	1
4		1098707	Hour Meter Assembly with Clip	1	1
5		1139950	Switch, Powe On/Off	1	1
6		1130062-NLA	NLA - Label, Control Panel - 10LX	1	1
6		1125073	Label, Control Panel - 10LXO2	1	1
7	D,E	1128739-NLA	NLA - Flowmeter without Palnuts	1	1
7	D	1135961	Flowmeter Assembly	1	1
8		1180107	Connector, Barbed, Brass (1/4")	1	1
9		1128738	Label, Logo	1	1
10		1102696	Screw, Phillips Pan Head SEMS (#8-32 x 3/8")	1	2
11		1085200	Transformer Assembly	1	1
12	A,H	1110779-NLA	NLA - Tubing Assembly, Check Valve Discharge	1	1
13	H	1108815-NLA	NLA - Tubing, Silicon (1/4 x 16-1/4")	1	1
14		1108809	Tubing, PVC (1/4 x 12-3/4")	1	1
15		1128737	PCB Control Assembly - 10LX	1	1
15	I	1194797	PCB Control Assembly, 10LXO2	1	1
15	J	1128742	PCB Control Assembly - 10LXO2	1	1
15		1128742E	WARR RPLMT - PCB Control Assembly - 10LXO2	1	1
16		1108737	Harness, Wire, 10 Pin	1	1
17		1076275	Tubing, Tygon (1/8 x 18")	1	1
18		2001902-NLA	NLA - Screw, Phillips Pan Head (#4-40 x 1/4")	1	2
19		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
19		2003065-NLA	NLA - Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	2
20	C	1128379	Package, Knob, Flowmeter, Silver w/ Black Swish (Pkg. of 10)	10	1
20	B	1128380	Package, Knob, Flowmeter, Black w/ Silver Swish	10	1
21		2001764	Jumper Wire Assembly	1	1
22	F,G	2001682	Tie Wrap (8" L)	1	2
23	F,G	1110473	Tubing, PVC (1/4 x 1-3/4")	1	2
24	F,G	2000436	Check Valve, Barb (5/16")	1	1
25	F,G	P0750X030	Filter, HEPA Outlet	1	1
26	F,G	1143144	Tubing, PVC (1/4 x 5")	1	1
27		1081276	Control Panel Only	1	1
28		1112562-NLA	NLA - Tech Only - Seal, Circuit Breaker	1	1
29		1194537	Washer, Foam	1	1


Control Panel and P.C. Board

Platinum 10

NOTE: A - Includes Tubing, Check Valve, Elbow and Tie Wrap.
B - Knob is black with silver "swish". Used on Flowmeters with a Slotted Head Screw on top.
C - Knob is silver with black "swish". Used on Flowmeters with an Indented Phillips Head Screw on top.
D - When replacing Flowmeter Assembly Part #1135961 in units built (Before 4/06), it is required that the Patient Outlet Pressure be increased to 9psi. Contact Invacare Technical Support at 1-800-832-4707 for more details.
E - Flowmeter 1128739 is no longer available, refer to Part #1135961 for replacement.
F - "After 7/19/06" the location of the HEPA Filter has changed. The HEPA Filter is now located inside the cabinet.
G - For use with Control Panel Assembly (After 7/19/06).
H - For use with Control Panel Assembly (Before 7/20/06).
I - For use with Control Panel Assembly built (After 9/30/16).
J - For use with Control Panel Assembly built (Before 10/01/16).

Control Panel and P.C. Board

Perfecto2, Perfecto2V, PerfectoW



SRE_020A


Control Panel and P.C. Board

Perfecto2, Perfecto2V, PerfectoW

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1143576-NLA	NLA - Control Panel Assembly - P2/PW	1	1
1		1143553-NLA	NLA - Control Panel with SensO2 Assembly - PO2/PO2W	1	1
1		1164802	Control Panel w/ SensO2 Assembly - PO2V	1	1
2		1110561	Flowmeter without Palnuts - 5 Liter	1	1
3		1180107	Connector, Barbed, Brass (1/4")	1	1
4		1194537	Washer, Foam	1	1
5		1143503	Horn/Alarm Assembly	1	1
6		1098707	Hour Meter Assembly with Clip	1	1
7		1143475	PCB Assembly - P2/PW	1	1
7		1143486	PCB Assembly - PO2/PO2W	1	1
7	C	1194768	PCB Assembly - PO2V	1	1
7		1143486E	WARR RPLMT - PCB Assembly - PO2/PO2W	1	1
8		1143499	Screw, Phillips Hex Washer Head (#8-16 x 3/8)	1	2
9		1143542	Breaker, Circuit 6 AMP	1	1
10		2001682	Tie Wrap (8" L)	1	2
11		1110473	Tubing, PVC (1/4 x 1-3/4")	1	2
12		2000436	Check Valve, Barb (5/16")	1	1
13		1143144	Tubing, PVC (1/4 x 5")	1	1
14		P0750X030	Filter, HEPA Outlet	1	1
15		2000507	Palnut (9/16-18)	1	3
16		1139950	Switch, Powe On/Off	1	1
17		1008547	Cord, Elastic	1	1
18		1006525	Tie Wrap (3-29/32" L)	1	1
19		1095497	Decal, Invacare Logo (1.6)	1	1
20		1143545	Label, LED - P2/PW	1	1
20		1143479	Label, LED - PO2/PO2V/PO2W	1	1
21		1143469	Panel, Control	1	1
22		2001764	Jumper Wire Assembly	1	1
23		1163429	Package, Label, No Smoking, Perfecto2/Perfecto2V/PerfectoW	20	1
24		1112562-NLA	NLA - Tech Only - Seal, Circuit Breaker	1	1
26		1128379	Package, Knob, Flowmeter, Silver w/ Black Swish (Pkg. of 10)	10	1
26		1128380	Package, Knob, Flowmeter, Black w/ Silver Swish	10	1
27		2000799	Adapter with Tubing for Humidifier Bottle	1	1
28		1143501-NLA	NLA - Power Jumper Assembly	1	1
	A,B	1143477	Main Harness Assembly	1	1
	A,B	1143526	Wire Harness Valve Assembly	1	1
NOTE: A - Not Shown. B - Connects to the PC Board. C - 1164808 replaced with 1194768 and can be used on all PO2V models.					

Sound Box Assembly

All Models Except Platinum 5



Sound Box Assembly

All Models Except Platinum 5

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1073613-NLA	NLA - Gasket, Foam	1	1
2		2001361	Kit, Exhaust Fan Service	1	1
			Fan, Exhaust	1	1
			NLA - Mount, Cable Tie	1	4
			Tie Wrap (3-29/32" L)	1	1
			NLA - Screw, Phillips Hex Washer Head (#7-19 x 3/8")	1	1
			NLA - Screw, Phillips Pan Head (10-24 x 3/4")	1	4
			Tie Wrap (15-1/4" L)	1	1
			Tie Wrap (8" L)	1	4
			NLA - Instruction Sheet	1	1
3		2000765-NLA	NLA - Spring, Extension (7/16 x 7-1/2" L)	1	1
4		2001375	Hose, Reinforced Silicon (1/2 x 8-1/2")	1	1
5		P0876X026	Clamp, Hose	1	2
6		2000755-NLA	NLA - Mount, Cable Tie	1	4
7		1076750-NLA	NLA - Screw, Phillips Pan Head (10-24 x 3/4")	1	4
8		2001682	Tie Wrap (8" L)	1	10
9	K	1085014-NLA	NLA - Exchanger, Heat	1	1
9		2000884-NLA	NLA - Kit, Heat Exchanger	1	1
9			NLA - Exchanger, Heat	1	1
5			Clamp, Hose	1	2
8			Tie Wrap (8" L)	1	4
10		2000411-NLA	NLA - Rod, Support	1	2
11		P0870X050-NLA	NLA - Nut, Keps (#6-32)	1	1
12		1088281-NLA	NLA - Tubing, Silicon (1/2 x 5")	1	1
13		2001194-NLA	Spring, Compression (17-1/4")	1	1
14		P0875X720	NLA - Nut (3/4-16)	1	1
15		2002997-NLA	NLA - Label, Diagnostic - IRC5LX02	1	1
16		2001096-NLA	NLA - Fitting, Bulk Head (1/2")	1	1
17	J	1108078-NLA	NLA - Filter, Inlet	1	1
17		1131249	Filter, Inlet Sound Reduced	1	1
18		2002264	Tubing, Reinforced Braided (1/2 x 26")	1	2
19	A,C	1105552	Kit, 4 Way Valve/Manifold Service	1	1
			NLA - Manifold Assembly	1	1
			Kit, 4-Way Valve w/ Screws	1	1
			Hose, Reinforced Silicon (1/2 x 8-1/2")	1	1
			Tubing, Reinforced Braided (1/2x 4-1/8")	1	2
			Clamp, Hose	1	1
			Tubing, Reinforced Braided (1/2 x 26")	1	2
			Tie Wrap (8" L)	1	6
			Clamp, Plastic (1/2")	1	2
			Bushing, Step (1/2")	1	4
			Washer, Flat (#8)	1	2

Sound Box Assembly

All Models Except Platinum 5

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
			NLA - Screw, Phillips Pan Head (#8-32 x 3/4")	1	2
			NLA - Resistor Service Assembly	1	1
			Tie Wrap, Banded (3-29/32" L)	1	2
			Instruction Sheet	1	1
20	D	1083005-NLA	NLA - Valve, 4-Way with Screen	1	1
20	D	1083005E-NLA	NLA - Valve, 4-Way with Screws (Reconditioned)	1	1
20	E	1101141	Kit, 4-Way Valve w/ Screws	1	1
	B	1101141R	Valve, Reconditioned	1	1
	B	1101141E-NLA	Kit, Valve - Warranty Replacement	1	1
			Valve, Reconditioned	1	1
			Gasket	1	1
21		1095209	Screw, Socket Head (M5 x .8 x 2mm)	1	4
22	F	1095538-NLA	Tubing, Reinforced Braided (1/2x 4-1/8")	1	2
			NLA - Kit, Muffler Retro-Fit	1	1
			NLA - Muffler Assembly	1	1
			NLA - Screw, Phillips Pan Head (#10-16 x 1/2")	1	2
			Screw, Socket Button Head (1/4-20 x 3/8")	1	1
			NLA - Bracket, Mounting	1	1
			Clip, Flange	1	1
			NLA - Instruction Sheet	1	1
22	G	1085872-NLA	NLA - Muffler Assembly	1	1
23		2001116	Uni-Muffler	1	1
24		2001765	Harness, Horn/Alarm	1	1
25	L	2003052	Valve, Operator Pilot	1	1
		1085112	O-Ring, Pilot Valve (35/64" I.D.)	1	1
	B	1080973-NLA	NLA - Kit, Pilot Valve Stem and Poppet Replacement	1	1
			NLA - Pilot Valve Poppet and Stem	1	1
			NLA - Instruction Sheet, Pilot Coil	1	1
	B,H	2003083	Kit, Pilot Valve Poppet Replacement	1	1
	B,I	2003082-NLA	NLA - Instruction Sheet, Pilot Valve Poppet Kit	1	1
	B	1080971-NLA	NLA - Kit, Pilot Coil Replacement	1	1
			NLA - Coil, Pilot	1	1
			NLA - Instruction Sheet, Pilot Coil	1	1
	B	2000830	Guard, Cooler	1	1


Sound Box Assembly

All Models Except Platinum 5

NOTE: A - Manifold #2002752, 4-Way Valve #1083005 and 4-Way Valve #1083005E are no longer available. Must order Manifold/Valve Service Kit (p/n 1105552). To be used on units built before 01/01.
B - Not Shown.
C - Includes items 4, 5, 8, 18, 20, 21 and Instruction Sheet.
D - Before 12/00.
E - After 11/00.
F - Before 7/00.
G - After 6/00.
H - Includes O-Ring, Labels, Slotted Poppets and Instruction Sheet.
I - See Tech Services for a copy of Instruction Sheet.
J - Item no longer available, refer to Part #1131249.
K - Part Number 1085014 is No Longer Available by itself. Refer to the Heat Exchanger Kit part number 2000884
L - A new Pilot Valve O-Ring part #1085112 is recommended when ordering a new Operator Pilot Valve

Sound Box Assembly

Platinum 5 Models



Sound Box Assembly Platinum 5 Models


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		2001361	Kit, Exhaust Fan Service	1	1
			Fan, Exhaust	1	1
			NLA - Mount, Cable Tie	1	4
			Tie Wrap (3-29/32" L)	1	1
			NLA - Screw, Phillips Hex Washer Head (#7-19 x 3/8")	1	1
			NLA - Screw, Phillips Pan Head (10-24 x 3/4")	1	4
			Tie Wrap (15-1/4" L)	1	1
			Tie Wrap (8" L)	1	4
			NLA - Instruction Sheet	1	1
2		1108783	Tubing, Braided (1/2 x 23")	1	1
3		1108784	Tubing, Braided (1/2 x 20-1/4")	1	1
4		2001765	Harness, Horn/Alarm	1	1
5	B	1108085	Manifold/Valve Assembly	1	1
6		2001682	Tie Wrap (8" L)	1	5
7		1110848-NLA	NLA - Tubing, Braided (1/2 x 6-3/16")	1	2
8	F	1085872-NLA	NLA - Muffler Assembly	1	1
9	F	2001116	Uni-Muffler	1	1
10	M	P0875X560-NLA	NLA - Boot, Terminal	1	1
11	A	1110776-NLA	NLA - Capacitor Mounting Hardware	1	1
12	H	1068698-NLA	NLA - Tech Only - Capacitor, 15 MFD	1	1
12	D,I	1131250	Kit, Capacitor (10 MFD)	1	1
			NLA - Tech Only - Capacitor, 10 MFD	1	1
			Tie Wrap (17-1/2" L)	1	1
			Tie Wrap (8" L)	1	2
			Instruction Sheet, Capacitor Replacement	1	1
13		P0876X026	Clamp, Hose	1	5
14		1108822	Elbow, PVC Formed (1/2 x 2-1/2 x 2-1/2")	1	1
15	C	1106698-NLA	NLA - Elbow, Silicon (1/2 x 3-1/2")	1	1
16	J	1108817	Heat Exchanger Assembly	1	1
17		1110780-NLA	NLA - Stand-Off, Heat Exchanger	1	1
18		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	7
19		1108086-NLA	NLA - Tubing (1/2 x 19")	1	1
19		1112561	Tube, Compressor	1	1
20		1101141	Kit, 4-Way Valve w/ Screws	1	1
	L	1101141R	Valve, Reconditioned	1	1
	L	1101141E-NLA	Kit, Valve - Warranty Replacement	1	1
			Valve, Reconditioned	1	1
			Gasket	1	1
			Screw, Socket Head (M5 x .8 x 2mm)	1	4
21		1110852	Screw, Slotted Hex Washer Head (#10-24 x 1/2")	1	4
22		1003363-NLA	NLA - Screw, Phillips Pan Head (#8-32 x 3/8")	1	1
22		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	1

Sound Box Assembly Platinum 5 Models

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
23		1106675	Coupling, Inlet Filter	1	1
24		2000812	Seal, Foam	1	2
25	K	P0750X030	Filter, HEPA Outlet	1	1
26		1110777-NLA	NLA - Partition	1	1
26	E	1115183	Partition Assembly	1	1
27		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	1
28		1108815-NLA	NLA - Tubing, Silicon (1/4 x 16-1/4")	1	1
29		1108078-NLA	NLA - Filter, Inlet	1	1
29		1131249	Filter, Inlet Sound Reduced	1	1
30		1110768-NLA	NLA - Tech Only - Adaptor, Humidifier Bottle	1	1
30		2000799	Adapter with Tubing for Humidifier Bottle	1	1
31	F	1087225-NLA	NLA - Screw, Phillips Pan Head (#10-16 x 1/2")	1	2
32	G	1110775	Muffler/Y-Tube Assembly	1	1
	L	2003052	Valve, Operator Pilot	1	1
NOTE: A - Includes Screw, Tie Wraps and Foam Tape. B - Due to style change effective 7/24/03, includes Items 2, 3, 6, 13, 14, 20 and 32. C - Item no longer available, use Item 16 (1108817) as its replacement. D - Includes: (1) Capacitor, 10MFD and (2) Tie Wraps. E - After 12/5/02. F - Before 7/24/03. G - After 7/23/03. H - Use with Thomas Compressors. I - Use with Gast Compressors. J - Includes Heat Exchanger, Standoff, and Tie Wrap. K - After 7/19/06 the Hepa Filter was moved to the inside of cabinet, located between the Flow Meter and the Patient Outlet Barb. See "Control Panel and P.C. Plat. 5 - After 12/5/02" category to view reflective changes. L - Not Shown. M - Boot removed from assembly after 6/28/07.					

Sound Box Assembly

Platinum 10



Sound Box Assembly

Platinum 10

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		2000817	Fan, Exhaust	1	1
2		1110852	Screw, Slotted Hex Washer Head (#10-24 x 1/2")	1	4
3		1108783	Tubing, Braided (1/2 x 23")	1	1
4		1108784	Tubing, Braided (1/2 x 20-1/4")	1	1
5		2001765	Harness, Horn/Alarm	1	1
6	B,C	1130671	Manifold/Valve Assembly - 10 Liter	1	1
	H	2003052	Valve, Operator Pilot	1	1
		1085112	O-Ring, Pilot Valve (35/64" I.D.)	1	1
7		1101141	Kit, 4-Way Valve w/ Screws	1	1
	D	1101141R	Valve, Reconditioned	1	1
	D	1101141E-NLA	Kit, Valve - Warranty Replacement	1	1
			Valve, Reconditioned	1	1
			Gasket	1	1
			Screw, Socket Head (M5 x .8 x 2mm)	1	4
8		2001682	Tie Wrap (8" L)	1	5
9		1130670	Muffler/Y-Tube Assembly	1	1
10	E	P0875X560-NLA	NLA - Boot, Terminal	1	1
11		1128740	Capacitor, 30 MFD	1	1
12		1089455	Tie Wrap (15-1/4" L)	1	2
13		1006525	Tie Wrap (3-29/32" L)	1	3
14		P0876X026	Clamp, Hose	1	2
15	A	1108817	Heat Exchanger Assembly	1	1
16		1110780-NLA	NLA - Stand-Off, Heat Exchanger	1	1
17	G	1107422-NLA	NLA - Exchanger, Heat	1	1
18		1108822	Elbow, PVC Formed (1/2 x 2-1/2 x 2-1/2")	1	1
19	G	1106698-NLA	NLA - Elbow, Silicon (1/2 x 3-1/2")	1	1
20		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	7
21		1112561	Tube, Compressor	1	1
22		1003363-NLA	NLA - Screw, Phillips Pan Head (#8-32 x 3/8")	1	1
23		1106675	Coupling, Inlet Filter	1	1
24		2000812	Seal, Foam	1	2
25	F	P0750X030	Filter, HEPA Outlet	1	1
26		1115183	Partition Assembly	1	1
27		2003065-NLA	NLA - Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	1
27		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	1
28		1108815-NLA	NLA - Tubing, Silicon (1/4 x 16-1/4")	1	1
29		1131249	Filter, Inlet Sound Reduced	1	1
30		2000799	Adapter with Tubing for Humidifier Bottle	1	1


Sound Box Assembly

Platinum 10

NOTE: A - Includes Heat Exchanger, Clamp, Elbow, Standoff and Tie Wrap.
B - Includes Items 3, 4, 7-9, 14, 18, Manifold and Tie Wraps.
C - Includes O-Ring and Operator Pilot Valve effective 6/4/05.
D - Not Shown.
E - Boot removed from units after 6/28/07.
F - After 7/19/06 the Hepa Filter was moved to the inside of cabinet, located between the Flow Meter and the Patient Outlet Barb. See "Control Panel and P.C. 10 Liter Models" category to view reflective changes.
G - Item no longer available, use Item 11 (1108817) as its replacement.
H - A new Pilot Valve O-Ring part #1085112 is recommended when ordering a new Operator Pilot Valve

Sound Box Assembly

Perfecto2, Perfecto2V




Sound Box Assembly Perfecto2, Perfecto2V

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1131249	Filter, Inlet Sound Reduced	1	1
2		1143449	Restrictor, Inlet Filter	1	1
3		1106675	Coupling, Inlet Filter	1	1
4		1143435	Resonator Assembly	1	1
5		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	7
6		1143581	Tube, Compressor Inlet	1	1
7		2001682	Tie Wrap (8" L)	1	3
8		1143484	Fan, Exhaust - P2/PW	1	1
9		1143574-NLA	NLA - Transformer Assembly - P2/PW	1	1
9		1143478	Transformer Assembly - PO2/PO2V/PO2W	1	1
10		1143491	Rivet, Isolator Mount	1	4
11		1143554	Muffler Assembly	1	1
12		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	5
13		1147917	Mount, Foam Valve 4-Way	1	1
14		1143467	Foam, Soundbox - Vertical	1	1
15		1143474	Gusset, Soundbox	1	2
16	A,B	1172007-NLA	NLA - Heat Exchanger Assembly	1	1
16		1186169	Heat Exchanger Assembly	1	1
17		1143582	Clamp, Hose Adjustable (3/8")	1	1
18	C	1169408-NLA	NLA - Fitting, Elbow Swivel (1/4 NPT x 3/8 Tube)	1	1
18		1143489	Elbow, Male Flare	1	1
NOTE: A - Use this part number to replace the Heat Exchanger after 05/17/11. When ordering a new Heat Exchanger, you must also order a new fitting (1169408). B - No longer available as a service part. Use 1143557 to replace 1172007, you must also order 1143489 also. C - No longer available. Use 1143489 as its replacement.					

Sound Box Assembly

PerfectoW



Sound Box Assembly

PerfectoW


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1131249	Filter, Inlet Sound Reduced	1	1
2		1153703-NLA	NLA - Tube, Clear (1/2 x 2-3/4")	1	1
3		1106675	Coupling, Inlet Filter	1	1
4		1143435	Resonator Assembly	1	1
5		1128404	Screw, Phillips Hex Washer Head (#8-32 x 3/8")	1	7
6		1143473	Soundbox, Horizontal	1	1
7	A	2001682	Tie Wrap (8" L)	1	3
8		1143484	Fan, Exhaust - P2/PW	1	1
9		1143574-NLA	NLA - Transformer Assembly - P2/PW	1	1
9		1143478	Transformer Assembly - PO2/PO2V/PO2W	1	1
10		1143491	Rivet, Isolator Mount	1	4
11		1143554	Muffler Assembly	1	1
12		1108810	Screw, Hi-Lo Six Lobe Drive (#10-16 x 5/8")	1	5
13		1147917	Mount, Foam Valve 4-Way	1	1
14		1154030	Foam, Soundbox - Vertical	1	1
15		1143474	Gusset, Soundbox	1	2
16	B,C	1172007-NLA	NLA - Heat Exchanger Assembly	1	1
16		1186169	Heat Exchanger Assembly	1	1
17		1143582	Clamp, Hose Adjustable (3/8")	1	1
18		1153693	Foam, Resonator - Top	1	1
19		1153694	Foam, Resonator - Bottom	1	1
20		1153699	Foam, Soundbox - Ceiling	1	1
21		1153707	Tube, Resonator Inlet Preformed	1	1
22		1153692-NLA	NLA - Tube, Inlet Silencer	1	1
23		1153788	Resonator, Inlet	1	1
24		1153708-NLA	NLA - Tube, Resonator Outlet Preformed	1	1
25		1143445-NLA	NLA - Soundbox, Vertical	1	1
26	D	1169408-NLA	NLA - Fitting, Elbow Swivel (1/4 NPT x 3/8 Tube)	1	1
26	B	1143489	Elbow, Male Flare	1	1

NOTE: A - After 1/31/10, this number is the replacement for the numbers 1006525 and 2001694.
 B - Use this part number to replace Heat Exchanger after 05/17/11. When ordering a new Heat Exchanger, you must also order the fitting (1169408).
 C - No longer available as a service part. Use 1143557 to replace 1172007, you must also order 1143489.
 D - No longer available. Use 1143489 as its replacement.


Sieve Bed Assembly and Product Tank

All Models Except Platinum 5


**Models 5LX, 5LX02, 5LX02T
6LX02, 6LX02T, 6LX02H**



Before 8/13/98




After 8/12/98



Side View


Models 6LX02H 6LX02HT

From Control Panel




View of Product Tank
Sens02 Units

Model 3LX



From Control Panel



View of Product Tank
Standard 3LX, 5LX

Sieve Bed Assembly and Product Tank

All Models Except Platinum 5

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,P	1101707	Sieve Bed Assembly (Pair)	2	1
2		2001682	Tie Wrap (8" L)	1	1
3		2000897	Tubing, Tygon (1/4 x 24")	1	1
4		1104310	Kit, Check Valve	1	1
			Clamp, Hose	1	2
			Tie Wrap (3-29/32" L)	1	1
			Instruction Sheet	1	1
			Check Valve, Barb (5/16")	1	2
			Tie Wrap (8" L)	1	9
			Tube, Clear (1/4")	1	1
5	L,N	1182123	Kit, Regulator and Adaptor, No-Bleed Regulator	1	1
			Regulator/Adaptor Assembly, No-Bleed Regulator	1	1
			Screw, Regulator (8-32 x 3/8")	1	4
			NLA - O-Ring (11/64" I.D.)	1	2
5	D,G E,M,O	1075498-NLA	NLA - Regulator	1	1
5		1182024		1	1
			Adaptor, Regulator	1	1
			Regulator (5L)	1	1
			Restrictor (1/8" Barb)	1	1
			Tubing, Silicone (1/8")	1	1
			Tubing, PVC (1/4")	1	1
			Screw, Regulator (8-32 x 3/8")	1	4
			NLA - O-Ring (11/64" I.D.)	1	2
			Tie Wrap (3-29/32" L)	1	4
6	F,K	1082588-NLA	NLA - Product Tank Assembly - IRC5LX	1	1
6		1114925-NLA	NLA - Product Tank Assembly - IRC5LX	1	1
6	G,I	1075665-NLA	NLA - Product Tank Assembly - IRC5LX02	1	1
7		1002807	Clamp, Hose	1	2
8	H B B B	2001656	PE Valve Assembly w/ Fitting	1	1
		1080972	Kit, PE Valve Stem and Poppet Replacement	1	1
			PE Valve Poppet and Stem	1	1
			Kit, PE Coil/Valve Instruction Sheet	1	1
		1080970	NLA - Kit, PE Coil Replacement	1	1
			Coil, PEI	1	1
			Kit, PE Coil/Valve Instruction Sheet	1	1
		1081208	Kit, PE Coil/Valve Instruction Sheet	1	1
9		2000492	Elbow, Barbed (1/2")	1	2
10		2001652	Clamp, Plastic (1/2")	1	2
11		2002264	Tubing, Reinforced Braided (1/2 x 26")	1	1
12	D	1070147-NLA	NLA - Screw, Phillips Pan Head (#4-24 x 1")	1	4
12	E	1075666	Screw, Regulator (8-32 x 3/8")	1	4
13		2001321-NLA	NLA - T-Fitting	1	1
14		1075499-NLA	NLA - Regulator, Two Way	1	1

Sieve Bed Assembly and Product Tank


All Models Except Platinum 5

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
15		1076916-NLA	NLA - Cap, End, Vinyl, Black	1	1
16		1070050	Tie Wrap (3-29/32" L)	1	2
17		2001951-NLA	NLA - Restrictor Assembly - IRC3LX	1	1
17	C	1078263-NLA	NLA - Restrictor Assembly (IRC5LX)	1	1
18		1089455	Tie Wrap (15-1/4" L)	1	4
	B,D	2002808-NLA	NLA - Plate, Adaptor	1	1
	B,J	2001673-NLA	NLA - Kit, Pressure Gauge Tap In with Gauge	1	1
			NLA - Fitting, Tee Barb (1/4")	1	1
			NLA - Plug, Hose (1/4")	1	1
			NLA - Clamp (1/4")	1	4
			NLA - Tubing, Silicone (1/4 x 16")	1	1
			NLA - Kit, Pressure Gauge	1	1
			Instruction Sheet	1	1
	B	2001671	NLA - Kit, Pressure Gauge Tap In without Gauge	1	1
			NLA - Fitting, Tee Barb (1/4")	1	1
			NLA - Plug, Hose (1/4")	1	1
			NLA - Clamp (1/4")	1	4
			NLA - Tubing, Silicone (1/4 x 16")	1	1
			Instruction Sheet	1	1
		2001674-NLA	NLA - Kit, Pressure Gauge	1	1

NOTE: A - Includes (2) Sieve Beds, (2) Clamps, (2) Hose Clamps and (4) Tie Wraps.
 B - Not Shown.
 C - Used on Serial Numbers: 99A655582, 99A655585, 99A655590, 99A655594, 99A655596, 99A65662, 99A65664, 99A65688, 99A65692, 99A65744, 99A65745, 99A65877, 99A65879, 99A65900.
 D - Before 5/98.
 E - After 4/98.
 F - LX Models.
 G - LX02 Models.
 H - Includes Tubing, Tie Wraps, Barb Fitting and PE Valve.
 I - Before 6/23/03.
 J - Item no longer available. Please order part #2001671 and 2001674 as a replacement.
 K - Item no longer available. Please order part #1114925 as a replacement.
 L - Kit includes the Regulator, (4) Screws, (2) O-Rings
 M - Kit includes Regulator, Adaptor, Restrictor, Screws, O-Rings, Tie Wraps, and Tubing.
 N - Used with IRC5LX.
 O - Used with IRC5LX02.
 P - May contain reactivated sieve.

Sieve Bed Assembly and Product Tank

Platinum 5 Models



Sieve Bed Assembly and Product Tank Platinum 5 Models

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,Q	1101707	Sieve Bed Assembly (Pair)	2	1
2		1110774	Tubing, PVC (1/4 x 2-1/2")	1	2
3		2001682	Tie Wrap (8" L)	1	9
4		2001656	PE Valve Assembly w/ Fitting	1	1
5		1108783	Tubing, Braided (1/2 x 23")	1	1
6		1108784	Tubing, Braided (1/2 x 20-1/4")	1	1
7		2001652	Clamp, Plastic (1/2")	1	2
8		1002807	Clamp, Hose	1	2
9		1104310	Kit, Check Valve	1	1
	B		Clamp, Hose	1	2
			Tie Wrap (3-29/32" L)	1	1
			Instruction Sheet	1	1
			Check Valve, Barb (5/16")	1	2
			Tie Wrap (8" L)	1	9
			Tube, Clear (1/4")	1	1
10		1182123	Kit, Regulator and Adaptor, No-Bleed Regulator	1	1
			Regulator/Adaptor Assembly, No-Bleed Regulator	1	1
			Screw, Regulator (8-32 x 3/8")	1	4
	L,O		NLA - O-Ring (11/64" I.D.)	1	2
10		1182500	Kit, Regulator/Adaptor w/ Bled Adaptor, Foam Mount and Check Valve Assy.	1	1
			Kit, Regulator and Adaptor w/ Bled Regulator	1	1
			Mount, Foam Product Tank	1	1
			Check Valve Assembly (14-1/2")	1	1
			Tie Wrap (17-1/2" L)	1	1
			Tie Wrap (8" L)	1	2
			Instruction Sheet, Contents and Application	1	1
			Adaptor, Regulator	1	1
	N,P		Regulator (5L)	1	1
			Restrictor (1/8" Barb)	1	1
			Tubing, Silicone (1/8")	1	1
			Tubing, PVC (1/4")	1	1
			Screw, Regulator (8-32 x 3/8")	1	4
			NLA - O-Ring (11/64" I.D.)	1	2
			Tie Wrap (3-29/32" L)	1	4
11		1116731-NLA	NLA - Product Tank Assembly - IRC5LX	1	1
11		1075665-NLA	NLA - Product Tank Assembly - IRC5LX02	1	1
11	E	1116730	Product Tank Assembly - IRC5LX02	1	1
12	D,F,G	1106774-NLA	NLA - HomeFill Tubing Assembly	1	1
12	F	1116433	HomeFill Tubing Assembly	1	1
12	C	1108809	Tubing, PVC (1/4 x 12-3/4")	1	1
13	J	1070050	Tie Wrap (3-29/32" L)	1	1
14					

Sieve Bed Assembly and Product Tank


Platinum 5 Models

NOTE: A - Includes (2) Sieve Beds, (2) Clamps, (2) Hose Clamps and (4) Tie Wraps.
B - Includes Check Valve, Tubing and Tie Wraps.
C - Includes Tubing and Fittings.
D - Includes Items 8 and 11.
E - Before 6/23/03.
F - After 6/22/03.
G - Individual item no longer available.
J - Due to design change the Homefill Tubing Assembly has changed from a 8-1/2" to a 18" Hose, with a removal of the fitting as well. The change is compatible with existing units.
L - Kit includes the Regulator, (4) Screws, (2) O-Rings
N - Kit includes Regulator, Adaptor, Tank Mount, Check Valve, Tie Wraps, and Instruction Sheet.
O - Used with IRC5LX.
P - Used with IRC5LXO2.
Q - May contain reactivated sieve.

Sieve Bed Assembly and Product Tank

Platinum 10

NOTE: Item #11 Product Tank is no longer available as an individual part. Please refer to item #10 Product Tank Assembly.



Sieve Bed Assembly and Product Tank

Platinum 10


NOTE: Item #11 Product Tank is no longer available as an individual part. Please refer to item #10 Product Tank Assembly.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	B	1144034	Kit, Sieve Beds (Pair)	1	1
			Sieve Bed Assembly - 10LX/10LXO2	1	2
			Tie Wrap (8" L)	1	4
			Clamp, Plastic (1/2")	1	2
			Clamp, Hose	1	2
2		1110774	Tubing, PVC (1/4 x 2-1/2")	1	2
3		2001682	Tie Wrap (8" L)	1	9
4	A	2001656	PE Valve Assembly w/ Fitting	1	1
5		1108783	Tubing, Braided (1/2 x 23")	1	1
6		2001652	Clamp, Plastic (1/2")	1	2
7		1002807	Clamp, Hose	1	2
8		1108784	Tubing, Braided (1/2 x 20-1/4")	1	1
9		1075666	Screw, Regulator (8-32 x 3/8")	1	4
10	C	1132389	Product Tank Assembly - 10LX	1	1
10	C	1130675	Product Tank Assembly - 10LXO2	1	1
12		1104310	Kit, Check Valve	1	1
			Clamp, Hose	1	2
			Tie Wrap (3-29/32" L)	1	1
			Instruction Sheet	1	1
			Check Valve, Barb (5/16")	1	2
			Tie Wrap (8" L)	1	9
			Tube, Clear (1/4")	1	1
13	H,K	1130067-NLA	NLA - Regulator, Manifold Mount, 10L w/o O2 Bleed	1	1
13	I,J	1130066-NLA	NLA - Regulator, Manifold Mount, 10L w/ O2 Bleed	1	1
14	E,G	1116433	HomeFill Tubing Assembly	1	1
15		1108809	Tubing, PVC (1/4 x 12-3/4")	1	1
16		1070050	Tie Wrap (3-29/32" L)	1	1

NOTE: A - Includes Valve, Tubing and Tie Wraps.
 B - Includes Item #1 (2) Sieve Bed Assemblies, item #3 (4) Tie Wraps, item #6 (2) Plastic Clamps and item #7 (2) Hose Clamps
 C - Includes Product Tank, Check Valve Assembly, Tie Wraps and Homefill Tubing. Item #11 Product Tank is no longer available as an individual part.
 E - Includes Fitting, Restrictor, Hose, Filter and Tie Wraps.
 G - Due to design change the Homefill Tubing Assembly has changed from a 8-1/2" to a 18" Hose, with a removal of the fitting as well. The change is compatible with existing units.
 H - Regulator does not have an O2 bleed port
 I - Regulator does have an O2 bleed port
 J - No longer available as a service part, please use 1130675 as its replacement.
 K - No longer available as a service part, please use 1132389 as its replacement

Sieve Bed Assembly and Product Tank

Perfecto2



SRE_002

Sieve Bed Assembly and Product Tank

Perfecto2

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1159071	Product Tank Assembly w/o SensO2 Regulator	1	1
1		1182228	Product Tank Assembly w/ SensO2 Regulator	1	1
2		2001682	Tie Wrap (8" L)	1	3
3	A	1182123	Kit, Regulator and Adaptor, No-Bleed Regulator	1	1
			Regulator/Adaptor Assembly, No-Bleed Regulator	1	1
			Screw, Regulator (8-32 x 3/8")	1	4
			NLA - O-Ring (11/64" I.D.)	1	2
3	B	1182500	Kit, Regulator/Adaptor w/ Bled Adaptor, Foam Mount and Check Valve Assy.	1	1
			Kit, Regulator and Adaptor w/ Bled Regulator	1	1
			Mount, Foam Product Tank	1	1
			Check Valve Assembly (14-1/2")	1	1
			Tie Wrap (17-1/2" L)	1	1
			Tie Wrap (8" L)	1	2
			Instruction Sheet, Contents and Application	1	1
4	C	1143544	Check Valve Assembly (9")	1	1
5	D	1182170	Check Valve Assembly (14-1/2")	1	1
6	E,I	1161257	Kit, Sieve Beds (Pair)	1	1
7			Tie Wrap (8" L)	1	2
8			Tie Wrap (17-1/2" L)	1	4
9			Sieve Bed Assembly	1	2
7		2001682	Tie Wrap (8" L)	1	2
8		1018239	Tie Wrap (17-1/2" L)	1	4
10		1143578	PE Valve Assembly	1	1
10		1149649	Kit, PE Valve w/ Fittings	1	1
			Valve, PE w/ Restrictor (24V Coil)	1	1
			NLA - Fitting, Barb (1/4 x 1/8")	1	2
11		1143447	Tube, Exhaust	1	1
12		2001682	Tie Wrap (8" L)	1	3
13		1143446	Tube, F	1	1
14	F	1143556-NLA	NLA - Tech Only - Manifold/Valve Assembly	1	1
	G	1101141	Kit, 4-Way Valve w/ Screws	1	1
	H	2003052	Valve, Operator Pilot	1	1
		1085112	O-Ring, Pilot Valve (35/64" I.D.)	1	1
15		P0876X026	Clamp, Hose	1	4
16		1188708	Tube, Polyurethane, Preformed (1/2")	1	2


Sieve Bed Assembly and Product Tank

Perfecto2

NOTE: A - Used with Product Tank Assembly (1159071), Kit includes the Regulator and Adaptor Assembly, (4) Screws, (2) O-Rings
B - Used with Product Tank Assembly (1182228), Kit includes Regulator and Adaptor Assembly, Foam Product Tank Mount, Check Valve Assembly, Instruction Sheet, and Tie Wraps.
C - Includes (1) Valve, (2) Tie Wraps, (1) 6" Tubing and (1) 3" Tubing.
D - Includes (1) Valve, (2) Tie Wraps, (1) 10" Tubing and (1) 4-1/2" Tubing.
E - Includes items 7, 8, and 9.
F - Includes (1) 4 Way Valve w/ Screws, (1) Pilot Valve O-Ring, (1) Air Manifold, (1) Operator Pilot Low Power Valve.
G - Includes the 4 Way Valve w/ Screws.
H - A new Pilot Valve O-Ring (1085112) is recommended when ordering 2003052.
I - May contain reactivated sieve.

Sieve Bed Assembly and Product Tank

Perfecto2V




Sieve Bed Assembly and Product Tank

Perfecto2V

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1182226	Product Tank Assembly w/ Senso2 Regulator	1	1
2		2001682	Tie Wrap (8" L)	1	3
3	A	1182500	Kit, Regulator/Adaptor w/ Bled Adaptor, Foam Mount and Check Valve Assy.	1	1
			Kit, Regulator and Adaptor w/ Bled Regulator Mount, Foam Product Tank	1	1
			Check Valve Assembly (14-1/2")	1	1
			Tie Wrap (17-1/2" L)	1	1
			Tie Wrap (8" L)	1	2
			Instruction Sheet, Contents and Application	1	1
4	B	1143544	Check Valve Assembly (9")	1	1
5	C	1182170	Check Valve Assembly (14-1/2")	1	1
6	D,I	1165099	Kit, Sieve Beds (Pair)	1	1
7			Tie Wrap (8" L)	1	4
8			Tie Wrap (17-1/2" L)	1	4
9			NLA - Tech Only - Sieve Bed Assembly	1	2
7		2001682	Tie Wrap (8" L)	1	4
8		1018239	Tie Wrap (17-1/2" L)	1	4
10	E	1165100	Kit, Adapters, Sieve Bed	1	1
11			NLA - Tech Only - Adapter, Sieve Bed Top	1	1
12			Adapter, Sieve Bed Bottom	1	1
13		1164785	PE Valve Assembly	1	1
14		1143447	Tube, Exhaust	1	1
15		2001682	Tie Wrap (8" L)	1	3
16		1143446	Tube, F	1	1
17	F	1143436	Manifold, Air	1	1
	G	1101141	Kit, 4-Way Valve w/ Screws	1	1
		2003052	Valve, Operator Pilot	1	1
		1085112	O-Ring, Pilot Valve (35/64" I.D.)	1	1
18	H	1160798-NLA	NLA - Clamp, Hose Spring (1/2")	1	4
18		P0876X026	Clamp, Hose	1	4
19		1188708	Tube, Polyurethane, Preformed (1/2")	1	2
NOTE: A - Kit includes the Regulator and Adaptor Assembly, Product Tank Foam Mount, Check Valve Assembly, Instruction Sheet, and Tie Wraps. B - Includes (1) Valve, (2) Tie Wraps, (1) 6" Tubing and (1) 3" Tubing. C - Includes (1) Valve, (2) Tie Wraps, (1) 10" Tubing and (1) 4-1/2" Tubing. D - Includes items 7-9 E - Includes items 11 and 12. Kit accommodates (1) Sieve Bed F - Order with 1101141, 2003052, and 1085112. G - Includes the 4-Way Valve with Screws only H - No Longer Available. Use part #P0876X026 as a replacement. I - May contain reactivated sieve.					

Sieve Bed Assembly and Product Tank

PerfectoW



Sieve Bed Assembly and Product Tank

PerfectoW

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1161712-NLA	NLA - Product Tank Assembly	1	1
1	B	1182226	Product Tank Assembly w/ Senso2 Regulator	1	1
2		2001682	Tie Wrap (8" L)	1	3
3	A,C	1182123	Kit, Regulator and Adaptor, No-Bleed Regulator	1	1
			Regulator/Adaptor Assembly, No-Bleed Regulator	1	1
			Screw, Regulator (8-32 x 3/8")	1	4
			NLA - O-Ring (11/64" I.D.)	1	2
3	B,D	1182500	Kit, Regulator/Adaptor w/ Bled Adaptor, Foam Mount and Check Valve Assy.	1	1
			Kit, Regulator and Adaptor w/ Bled Regulator	1	1
			Mount, Foam Product Tank	1	1
			Check Valve Assembly (14-1/2")	1	1
			Tie Wrap (17-1/2" L)	1	1
			Tie Wrap (8" L)	1	2
			Instruction Sheet, Contents and Application	1	1
4	E	1143544	Check Valve Assembly (9")	1	1
5	F	1182170	Check Valve Assembly (14-1/2")	1	1
6		1006525	Tie Wrap (3-29/32" L)	1	3
7	G,M	1161257	Kit, Sieve Beds (Pair)	1	1
8			Tie Wrap (8" L)	1	2
9			Tie Wrap (17-1/2" L)	1	4
10			Sieve Bed Assembly	1	2
8		2001682	Tie Wrap (8" L)	1	2
9		1018239	Tie Wrap (17-1/2" L)	1	4
11		1143578	PE Valve Assembly	1	1
11		1149649	Kit, PE Valve w/ Fittings	1	1
			Valve, PE w/ Restrictor (24V Coil)	1	1
			NLA - Fitting, Barb (1/4 x 1/8")	1	2
12		1143447	Tube, Exhaust	1	1
13		2001682	Tie Wrap (8" L)	1	3
14		1143446	Tube, F	1	1
15	H	1143556-NLA	NLA - Tech Only - Manifold/Valve Assembly	1	1
	I	1101141	Kit, 4-Way Valve w/ Screws	1	1
	J	2003052	Valve, Operator Pilot	1	1
		1085112	O-Ring, Pilot Valve (35/64" I.D.)	1	1
16	K	1160798-NLA	NLA - Clamp, Hose Spring (1/2")	1	4
16		P0876X026	Clamp, Hose	1	4
17		1188708	Tube, Polyurethane, Preformed (1/2")	1	2

Sieve Bed Assembly and Product Tank




PerfectoW

NOTE: A - Use with PW Models.
B - Use with PO2W Models.
C - Kit includes the Regulator/Adaptor Assembly, Screws, and O-Rings.
D - Kit includes Regulator/Adaptor Assembly, Product Tank Foam Mounting, Check Valve Assembly, Instruction Sheet, and Tie Wraps.
E - Includes Valve, (2) Tie Wraps, 6" Tubing and 3" Tubing.
F - Includes Valve, (2) Tie Wraps, 10" Tubing and 4-1/2" Tubing.
G - Includes items 8 - 10.
H - Includes 4 Way Valve w/ Screws, Pilot Valve O-Ring, Air Manifold, and Operator Pilot Low Power Valve.
I - Includes 4 Way Valve w/ Screws only.
J - A new Pilot Valve O-Ring (1085112) is recommended when ordering 2003052.
K - No Longer Available. Use part #P0876X026 as a replacement.
M - May contain reactivated sieve.


Remote Flowmeters

IRCRF5 - Adult, IRCPF16 - Pediatric

Model IRCRF5 - Adult



Model IRCPF16 - Pediatric



Remote Flowmeters

IRCRF5 - Adult, IRCPF16 - Pediatric

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		2001208	Handle, Flowmeter Stand, Chrome	1	1
2		2001206-NLA	NLA - Label, No Smoking	1	1
3	D	2003080	Stand, Flowmeter	1	1
4	B	2001207-NLA	NLA - Stand, Flowmeter	1	1
5		2001202-NLA	NLA - Fitting, Tee	1	1
6		1141758-NLA	NLA - Tubing, Silicone (1/4 x 16")	1	1
7		1073266-NLA	NLA - Hose, Silicone (1/8 x 6")	1	1
8		1180032	Hose, Silicone (1/4 x 1")	1	1
9		1073268-NLA	NLA - Hose, Silicone (1/8 x 2-13/16")	1	1
10		1073269-NLA	NLA - Hose, Silicone (1/8 x 3-5/16")	1	1
11	D	2003074	Mat, Rubber	1	1
11	B	1079652	Mat, Rubber	1	1
12		1026031-NLA	NLA - Nut, Hex with Nylon Cap (#10-32)	1	1
13		2001198-NLA	NLA - Fitting, Barbed (9/16-18 x 1/8")	1	1
14	D	2001103	Flowmeter	1	1
14	B,C	2001896-NLA	NLA - Flowmeter	1	1
15	D	2001288	Fitting, Barbed with Palnut (1/8")	1	1
15	B	1180107	Connector, Barbed, Brass (1/4")	1	1
16	D	2003075	Decal, Logo, Ped	1	1
16	B	2001210-NLA	NLA - Decal, Logo, In.	1	1
17		1018138	Screw, Socket Button Head (#10-24 x 1/4")	1	1
18		2000507	Palnut (9/16-18)	1	4
19		2001879-NLA	NLA - Restrictor, Barbed (1/8")	1	1
20		2001098	Tee, Single Barbed (1/8")	1	1
21		P0955X030-NLA	NLA - Connector, Humidifier Bottle	1	1
22		1067973	Screw, Phillips Flat Head (4-40 x 3/8")	1	4
23		3260S	Bottle, Humidifier	1	1
24		2000799	Adapter with Tubing for Humidifier Bottle	1	1
25		2000503	Adapter, Humidifier Bottle	1	1
	A	2001219	Carton	1	1
	A,D	2001105	Instruction Manual, Pediatric Flowmeter IRCPF16	1	1
	A,B	2001222-NLA	NLA - Instruction Sheet, Remote Flowmeter IRCRF5	1	1
	A,D	1071207	Tubing, Oxygen Supply (7' L)	1	1
	A,D	2003110	Tubing, Oxygen Supply (21' L)	1	1
	A,B	2001221-NLA	NLA - Tubing, Oxygen Supply (50' L)	1	1
	A,D	1071275-NLA	NLA - Kit, Set-Up Pediatric Flowmeter	1	1
			Tubing, Oxygen Supply (21' L)	1	1
			Tubing, Oxygen Supply (7' L)	1	2
			NLA - Connector, Tubing	1	1
			Adapter, Humidifier Bottle	1	1
			Instruction Manual, Pediatric Flowmeter IRCPF16	1	1

Remote Flowmeters


IRCRF5 - Adult, IRCPF16 - Pediatric

NOTE: A - Not Shown.
B - Adult.
C - Includes item 18.
D - Pediatric.


Invacare Check 02 Plus Concentrator Analyzer

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
		2002246-NLA	NLA - Tubing (12" L)	1	1


HomeFill Base and HomeFill Stand Components




Use on units with
straight tubing.



Use on units with
Coiled Tubing.



HomeFill Stand
Shown for
Reference Only



HomeFill Base and HomeFill Stand Components

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1099820	Fuse (2 AMP)	1	1
2		2000489	Filter, Cabinet	1	1
3	C	1108211	Cap, End, Dust, Green	1	1
4	F	1127653-NLA	NLA - Dust Cover with Lanyard (1-1/2")	1	1
4	M,N	1156644-NLA	NLA - Dust Cover with Lanyard (1")	1	1
4	M	60103657	Kit, Cap, Connector Dust, HomeFill II (Pkg. of 10)	10	1
5	B,E	1086773-NLA	NLA - Tip, Rubber, Black (3/4" Tubing)	1	4
5	E	1128290-NLA	NLA - Package, Rubber Tip, Black (Fits 3/4" Tubing)	4	1
	A,D	1157640	Package, Fitting Cover	1	1
	A	1112949	Cord, Power	1	1
6	K	1006525	Tie Wrap (3-29/32" L)	1	2
7	K	1145488	Package, Straight Internal Valve Fitting (1/8")	10	1
8	G,I	1105177	Supply Line Assembly w/ Coiled Tubing (12 ft.)	1	1
	A,H	1072004	Pad, Non-Skid Adhesive Back	1	1
9	I,L	1145487	Supply Line Assembly w/ Straight Tubing (3-21/64 FT)	1	1
10	J	1145489	Package, Hose Fitting (1/4" Barb)	10	1
11	J	2001682	Tie Wrap (8" L)	1	2

NOTE: A - Not Shown.
 B - Item no longer available, see Part #1128290.
 C - Used for Homefill Interface Port.
 D - For use with Homefill Oxygen Cylinder.
 E - For use on IOH250 Invacare Homefill Stand.
 F - For use over Oxy Lock Coupler with 1-1/2" Diameter.
 G - Includes Coiled Tubing, Fittings and Tie Wraps.
 H - Use with Homefill Ready Rack.
 I - After 4/12/06, the tubing style has changed to a Coiled Tubing.
 J - Use with Supply Line Assembly "Before 4/13/06" with Straight Tubing.
 K - Use with Supply Line Assembly "After 4/12/06" with Coiled Tubing.
 L - Includes Straight Tubing, Fittings and Tie Wraps.
 M - For use over Oxy Lock Coupler with 1" Diameter.
 N - See item #4 60103657 as a replacement for part# 1156644.

NOTE: Please contact an Invacare Parts Specialist for replacement parts at 1-800-333-6900 ext. 3722.

The diagram shows a green gas cylinder with a black pulse sensor attached to its top. The sensor is labeled 'Pulse Sensor' and 'Pulse Sensor'. The sensor is connected to a black hose. The diagram includes numbered callouts 1 through 6, which correspond to the following parts:





- 1: The black pulse sensor.
- 2: The black hose.
- 3: A close-up view of the sensor's connection point.
- 4: A close-up view of the sensor's connection point.
- 5: The gas cylinder.
- 6: The black hose.

**Cylinder
Shown for
Reference Only**

HomeFill Cylinder Assembly

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1157638	Kit, Protective Cap w/Hardware, Black	1	1
			Cap, End, Protective, Black	1	1
			Screw, Cap Mounting	1	3
2		1157640	Package, Fitting Cover	10	1
3	A	1117149	Kit, Receptacle Fitting Brass	1	1
			NLA - Fitting, Brass (1/4")	1	1
			Instruction Sheet	1	1
4	A	1157639	Kit, Fitting w/ Instructions	1	1
			Instruction Sheet	1	1
			Fitting, Male	1	1
5		1125090	Label, Oxygen Cylinder	1	1
6		1157637	Gauge, Pressure (1-1/4" Diam.)	1	1
NOTE: A - Existing Fitting shape must be determined before ordering replacement part, refer to image as a reference.					


Venture Demand Oxygen Delivery Device & Cylinder Accessories



Venture Demand Oxygen Delivery Device & Cylinder Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		2002246-NLA	NLA - Tubing (12" L)	1	1
2		2002919-NLA	NLA - Pouch, Venture	1	1
3	A	1147509	Bag Assembly, M2 Cylinder	1	1
4	B	1093320	Bag, Shoulder - M6/M9	1	1
5	C	1130794-NLA	NLA - Strap, Shoulder - ML6 and ML9 Shoulder Bag	1	1
5		1154103	Strap, Shoulder Bag	1	1
NOTE: A - Includes Bag and Strap. B - Bag color has changed from Green to Black, as of 02/10/06. C - No Longer Available. Replaced with 1154103.					


Concentrator Tool Kit



Concentrator Tool Kit

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1128756-NLA	Kit, Concentrator Tools	1	1
			NLA - Light, Pen-Size with Invacare Logo	1	1
			Screwdriver, Pen Size Phillips/Flat Head	1	1
			Tie Wrap (8" L)	1	2
			NLA - Multimeter with Leads	1	1
			Analyzer, Check O2 Plus Concentrator Invacare	1	1
			NLA - Strap, Ground with Wire and Clip	1	1
			NLA - Fitting, Tee Barb (1/4")	1	1
			NLA - Plug, Hose (1/4")	1	1
			NLA - Clamp (1/4")	1	4
			NLA - Tubing, Silicone (1/4 x 16")	1	1
			NLA - Kit, Pressure Gauge	1	1
			Instruction Sheet	1	1

Miscellaneous Components



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1168406	Kit, Knob, 0-6 lpm, Homefill Knob, 0-6 lpm, Homefill	1 10	1 1

Warranty Replacement Parts

Part Number	Description	Quantity per Assembly
1108720E-NLA See Note(s) A	WARR RPLMT - LXN PCB Assembly	1
1108731E-NLA See Note(s) A	WARR RPLMT - ISBN PCB Assembly	1
1106690E-NLA See Note(s) B	WARR RPLMT - Compressor Assembly, Thomas (5LX Plat.)	1
1128741E-NLA See Note(s) C	WARR RPLMT - Compressor (10LX/10LXO2)	1
1130666E See Note(s) G, C	WARR RPLMT - Compressor Assembly (10LX/10LXO2)	1
1143486E See Note(s) D	WARR RPLMT - PCB Assembly (PO2/PO2W)	1
1128742E See Note(s) E	WARR RPLMT - PCB Control Assembly (10LXO2)	1
1075071E-NLA See Note(s) F	WARR RPLMT - Compressor Assembly (6LX/6LXO2)	1
<p>NOTE: A - See "Control Panel and PC Boards (Platinum 5 Models - After 12/5/02)" Category for more information. B - See "Compressor and Mtg. Hardware (Platinum 5 Models)" Category for more information. C - See "Compressor and Mtg. Hardware (10 Liter Models)" Category for more information. D - See "Control Panel and PC Boards (Perfecto2/PerfectoW)" Category for more information. E - See "Control Panel - Rear View" Category for more information. F - See "Compressor and Mtg. Hardware (Non Platinum 5 Models)" Category for more information. G - When ordering 1130666E, you must order Kit part number 1186230.</p> <p>When ordering 1130666E, you must order Kit part number 1186230.</p>		

Invacare Corporation

www.invacare.com



One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Florida Technical Support
407-321-5630

Customer Service
800-333-6900

© 2020 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 03/3/20)

Form No. 93-38